

LW8372S

Sea Star Oocyte Immature

Whole immature sea star oocyte lysate boiled in SDS-PAGE sample buffer



KINEXUS

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Production

Product Name Long:	Sea star (<i>Pisaster ochraceus</i>) oocyte - immature - Whole lysate
Production Method:	Oocytes teased from ripe seastar oocyte gonads were resolved from follicle cells by 3 rounds of centrifugation at 2500 rpm for 5 min each time. Lysates were prepared from isolated oocytes that were homogenized by sonication in buffer formulated with 1% Triton X-100, 60 mM β -glycerophosphate, pH 7.2, 20 mM MOPS, 20 mM sodium pyrophosphate, 30 mM sodium fluoride, 5 mM EDTA, 3 mM benzamidine, 2 mM EGTA, 1 mM sodium orthovanadate, 1 mM phenylmethylsulfonylfluoride, 1 mM dithiothreitol, 10 μ M leupeptin, and 5 μ M pepstatin A. Lysates were further diluted in SDS-PAGE sample buffer at a final concentration of 2 mg/ml.
Amount:	200 μ g
Protein Concentration:	2 mg/ml
Storage Stability:	1 year at -20°C

Applications

Lysate Use Description:	For testing antibodies by immunoblotting.
--------------------------------	---

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-888-546-3987