

AB-NK008-4P

AurKA Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-serine/threonine kinase AurKA (Aurora A, AIK, STK15)



KINEXUS

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Target Protein

Name Long:	Aurora Kinase A (serine/threonine protein kinase 6)
Alias:	AIK; ARK1; AURKA; Aurora kinase A; Aurora,IPL1-related kinase 1; Aurora/IPL1-related kinase 1; Aurora-1; Aurora-A; Aurora-family kinase 1; STK7; BTAK; STK15; STK6; AURORA2; AurA; MGC34538; CCDS13451.1; ENSG00000087586
UniProt ID:	O14965
Sequence Predicted Mass (KDa):	45.809 (403 AA; O14965)
Observed SDS-PAGE Mass (KDa):	45-50

Immunogen

Antibody Immunogen Source:	Human AurKA (Aurora A, AIK, STK15) sequence peptide Cat. No.: PE-01ASI95
Antibody Immunogen Sequence:	CESKKRQWALED FEI
Location in Target:	Corresponds to amino acid residues E122 to I135
Peptide Type:	For pan-specific recognition of target expression levels.
Target Phosphosite:	Not phosphorylated

Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
Antibody Amount:	25 µg
Antibody Concentration:	1.0 mg/ml
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

AB-NK008-4P

AurKA Antibody



KINEXUS

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Applications

Product Use:	Western blotting Antibody microarrays
Antibody Dilution Recommended:	1 µg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Antibody Positive Controls:	Strong immunoreactivity with recombinant human AurKA on protein dot blots.
Antibody Cross Reactivities:	No immunoreactivity on protein dot blots with recombinant human AurKB and AurKC.

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-866-KINEXUS(546-3987)