# AB-NK020-1 CaMK2g-PAD Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-serine/threonine kinase CaMK2g

KINEXUS

Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

### **Target Protein**

Name Long:	Calcium/calmodulin-dependent protein-serine kinase 2 gamma
Alias:	Calcium/calmodulin-dependent protein kinase II gamma; CaM kinase II gamma subunit; CaMK2-gamma; CAMKG; CaMK-II gamma; CaMK-II gamma subunit; CaM-kinase II gamma; CaM-kinase II gamma chain; EC 2.7.11.17; Kinase CaMK2-gamma
UniProt ID:	Q13555
Sequence Predicted Mass (KDa):	65.242 (588 AA; Q13555-6); 62.609 (558 AA; Q13555); 62.165 (556 AA; Q13555-8); 61.368 (547 AA; Q13555-2); 60.743 (542 AA; Q13555-9); 59.760 (533 AA; Q13555-7); 59.607 (529 AA; Q13555-4); 59.038 (527 AA; Q13555-5); 58.963 (524 AA; Q13555-3); 55.961 (495 AA; Q13555-10)
Observed SDS-PAGE Mass (KDa):	55-60

#### Immunogen

Antibody Immunogen Source:	Human CaMK2g sequence peptide Cat. No.: PE-01BCJ90
Antibody Immunogen Sequence:	CDGIKGSTESSNTTTEDEDL
Location in Target:	Corresponds to amino acid residues D376 to L394; Before CaMAD, highly conserved in all CaMK2 isoforms
Peptide Type:	For pan-specific recognition of target expression levels.
Target Phosphosite:	Not phosphorylated

### Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
Antibody Amount:	25 μg
Antibody Concentration:	0.5 / 0.28 mg/ml
Lot Number:	10320 / 160804
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

# AB-NK020-1 CaMK2g-PAD Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

### Applications

Product Use:	Western blotting   Antibody microarrays
<b>Antibody Dilution Recommended:</b>	2 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals

This product is for in vitro research use only and is not intended for use in humans or animals.