# AB-NK021-1 CaMPK4-NT (CaMK4-2) Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-serine/threonine kinase CaMK4 (CaMPK4)

CaMK4 (CaMPK4)

Address: 8755 Ash Street, Suite 1

Vancouver, British Columbia, Email: info@kinexus.ca Canada V6P 6T3 Phone: 604-323-2547

### **Target Protein**

Name Long:	Calcium/calmodulin-dependent protein-serine kinase 4
Alias:	Brain Ca++-calmodulin-dependent protein kinase type IV; Calcium/calmodulin-dependent protein kinase IV; Calcium/calmodulin-dependent protein kinase type IV catalytic chain; Calspermin; CAM kinase- GR; CAM kinase IV; CAM kinase-GR; CaMK IV; CAMK4; Kinase CaMK4; CaMK-GR; CaMKIV; KCC4 CaMK IV; MGC36771; CCDS4103.1; ENSG00000152495
UniProt ID:	Q16566
Sequence Predicted Mass (KDa):	51.926 (473 AA; Q16566)
Observed SDS-PAGE Mass (KDa):	55-60

#### Immunogen

Antibody Immunogen Source:	Human CaMK4 (CaMPK4) sequence peptide
Antibody Immunogen Sequence:	CSSVTSSTENLVPDYWID
Location in Target:	Corresponds to amino acid residues C14 to D31; N-terminus
Peptide Type:	For pan-specific recognition of target expression levels.
Target Phosphosite:	Not phosphorylated

### Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced on an A.B.I. automated peptide synthesizer. The peptide was cleaved from the synthesis resin with HF and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by amino acid analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
Antibody Amount:	25 μg
Antibody Concentration:	NK021 (CaMK4-2) Lot 010320 with 0.5 mg/ml conc. ? mg/ml
Lot Number:	NK021 (CaMK4-2) Lot 010320 with 0.5 mg/ml conc. ?
Storage Buffer:	100 mM Tris-glycine, pH 7.0
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months, but either 0.1% sodium azide or 0.05% Thimerasol should be added. Avoid repeated freeze-thaw cycles.

# AB-NK021-1 CaMPK4-NT (CaMK4-2) Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

## Applications

•	
Product Use:	Western blotting   Antibody microarrays
<b>Antibody Dilution Recommended:</b>	2 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Antibody Positive Controls:	Strong immunoreactivity on Western blots with recombinant mouse calmodulin- dependent protein kinase 4 expressed in sf9 insect cells.

This product is for in vitro research use only and is not intended for use in humans or animals.