AB-NK052-4P EGFR Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-tyrosine kinase EGFR (ErbB1)



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Target Protein

Name Long:	Epidermal growth factor receptor-tyrosine kinase
Alias:	EGFR; Epidermal growth factor receptor; ErbB-1; ErbB, mENA; HER1; Receptor tyrosine-protein kinase ErbB-1; V-erb-b oncogene homologue; PIG61; CCDS5514.1; ENSG00000146648
UniProt ID:	P00533
Sequence Predicted Mass (KDa):	134.277 (1210 AA; P00533); 77.312 (705 AA; P00533-3); 69.228 (628 AA; P00533-4); 44.664 (405 AA; P00533-2)
Observed SDS-PAGE Mass (KDa):	160-180

Immunogen

Antibody Immunogen Source:	Human EGFR (ErbB1) sequence peptide Cat. No.: PE-01ATU90
Antibody Immunogen Sequence:	CEYLRVAPQSSEFIGA
Location in Target:	Corresponds to amino acid residues E1196 to A1210; C-terminus
Peptide Type:	For pan-specific recognition of target expression levels.
Target Phosphosite:	Not phosphorylated

Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
Antibody Amount:	25 μg
Antibody Concentration:	1 mg/ml
Lot Number:	130906
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

AB-NK052-4P EGFR Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications

Product Use:	Western blotting Antibody microarrays
Antibody Dilution Recommended:	2 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Antibody Positive Controls:	Medium immunoreactivity with recombinant human EGFR on protein dot blots.
Overall Antibody Specificity:	Very high selectivity
Antibody Cross Reactivities:	No immunoreactivity on protein dot blots with recombinant human ErbB2 and ErbB4.

This product is for in vitro research use only and is not intended for use in humans or animals.