

# AB-NK062-3

## FGFR1-3 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-tyrosine kinase FGFR1



**KINEXUS**

Address: 8755 Ash Street, Suite 1  
Vancouver, British Columbia,  
Canada V6P 6T3

Email: [info@kinexus.ca](mailto:info@kinexus.ca)  
Phone: 604-323-2547

### Target Protein

<b>Name Long:</b>	Fibroblast growth factor receptor-tyrosine kinase 1; Basic fibroblast growth factor receptor 1
<b>Alias:</b>	BFGFR; BFGF-R; CD331; CEK; C-fgr; FGFR; FGFR-1; Fgr; FLT2; KAL2; Pfeiffer syndrome; HBGR; KAL2; N-SAM; OGD; bFGF-R-1; H2; H3; H4; H5; ENSG00000077782
<b>UniProt ID:</b>	P11362
<b>Sequence Predicted Mass (KDa):</b>	95.344 (853 AA; P11362-21); 91.868 (822 AA; P11362); 91.760 (822 AA; P11362-19); 91.668 (820 AA; P11362-4); 91.580 (820 AA; P11362-14); 90.618 (812 AA; P11362-20); 82.162 (733 AA; P11362-6); 81.962 (731 AA; P11362-8); 81.875 (731 AA; P11362-15); 74.133 (662 AA; P11362-10); 73.933 (660 AA; P11362-12); 73.475 (662 AA; P11362-2); 73.274 (660 AA; P11362-5); 63.769 (573 AA; P11362-7); 63.569 (571 AA; P11362-9); 55.740 (502 AA; P11362-11); 55.540 (500 AA; P11362-13); 33.412 (302 AA; P11362-17); 33.125 (300 AA; P11362-18); 16.487 (150 AA; P11362-16); 6.682 (61 AA; P11362-3)
<b>Observed SDS-PAGE Mass (KDa):</b>	100-110

### Immunogen

<b>Antibody Immunogen Source:</b>	Human FGFR1 sequence peptide Cat. No.: PE-01AUD99
<b>Antibody Immunogen Sequence:</b>	CRHPAQLANGGLKRR
<b>Location in Target:</b>	Corresponds to amino acid residues R809 to R822; C-terminus
<b>Peptide Type:</b>	For pan-specific recognition of target expression levels.
<b>Target Phosphosite:</b>	Not phosphorylated

### Production

<b>Antibody Host Species:</b>	Rabbit
<b>Antibody Type:</b>	Polyclonal
<b>Antibody Ig Isotype Clone Lot:</b>	Immunoglobulin G
<b>Production Method:</b>	The immunizing peptide was produced by solid phase synthesis on a multi-peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
<b>Antibody Amount:</b>	25 µg
<b>Antibody Concentration:</b>	1 mg/ml
<b>Lot Number:</b>	140724
<b>Storage Buffer:</b>	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
<b>Storage Conditions and Stability:</b>	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

# AB-NK062-3

## FGFR1-3 Antibody



# KINEXUS

Address: 8755 Ash Street, Suite 1  
Vancouver, British Columbia,  
Canada V6P 6T3

Email: [info@kinexus.ca](mailto:info@kinexus.ca)  
Phone: 604-323-2547

### Applications

<b>Product Use:</b>	Western blotting   Antibody microarrays
<b>Antibody Dilution Recommended:</b>	2 µg/ml for immunoblotting
<b>Antibody Species Reactivity:</b>	Human, mouse, rat and many other mammals
<b>Antibody Positive Controls:</b>	Strong immunoreactivity with recombinant human FGFR1 on protein dot blots.
<b>Overall Antibody Specificity:</b>	High selectivity
<b>Antibody Cross Reactivities:</b>	No immunoreactivity on protein dot blots with recombinant human FGFR2, FGFR3 and FGFR4.

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit [www.kinexusproducts.ca](http://www.kinexusproducts.ca) or contact us at 1-866-KINEXUS(546-3987)