AB-NK065-2 Fyn-ANT Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-tyrosine kinase Fyn



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

	- 1	—	1 - " -
laro	ΔT	Lrα	TAIN
Targ			

Name Long:	Fyn proto-oncogene-encoded protein-tyrosine kinase
Alias:	FYN; FYN oncogene related to SRC, FGR, YES; Kinase Fyn; MGC45350; P59-FYN; P59-Fyn; SLK; SYN; hCG_34806; RP1-66H14.1-003; CCDS5094.1; ENSG00000010810
UniProt ID:	P06241
Sequence Predicted Mass (KDa):	60.762 (537 AA; P06241); 60.141 (534 AA; P06241-2); 54.513 (482 AA; P06241-3)
Observed SDS-PAGE Mass (KDa):	50-60

Immunogen

Antibody Immunogen Source:	Human Fyn sequence peptide Cat. No.: PE-01BCZ99
Antibody Immunogen Sequence:	CKDKEATKLTEERDGSLNQS
Location in Target:	Corresponds to amino acid residues C6 to S25; Near N-terminus
Peptide Type:	For pan-specific recognition of target expression levels.
Target Phosphosite:	Not phosphorylated

Production

Antibody Host Species:	Rabbit	
Antibody Type:	Polyclonal	
Antibody Ig Isotype Clone Lot:	Immunoglobulin G	
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.	
Antibody Amount:	25 μg	
Antibody Concentration:	0.5 mg/ml	
Lot Number:	160807	
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol	
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.	

AB-NK065-2 Fyn-ANT Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications

Product Use:	Western blotting Antibody microarrays
Antibody Dilution Recommended:	2 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals

This product is for in vitro research use only and is not intended for use in humans or animals.