

AB-NK073-3

ICK-1 Antibody



KINEXUS

Pan-specific rabbit polyclonal antibody for monitoring the expression of human protein-serine/threonine kinase ICK

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Target Protein

Protein Name Long:	Intestinal cell (MAK-like) kinase
Protein Alias:	ECO; HICK; Intestinal cell (MAK-like) kinase; Intestinal cell kinase; KIAA0936; Kinase ICK; Laryngeal cancer kinase 2; LCK2; MAK-related kinase; MGC46090; MRK; Serine/threonine kinase ICK
UniProt ID:	Q9UPZ9
Protein Molecular Mass:	71,427 Da (632 AA; Q9UPZ9); 34,071 Da (292 AA; Q9UPZ9-2)

Immunogen

Antibody Immunogen Source:	
Antibody Immunogen Sequence:	CGHPLGSTTQNLQDSEKPQK
Antibody Immunogen Description:	Corresponds to amino acid residues G287 to K305.
Antibody Target Type:	Pan-specific

Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Isotype:	IgG
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from these animals was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
Amount:	25 µg
Antibody Concentration:	0.9 mg/ml
Storage Buffer:	Phosphate buffered saline pH 7.4, 0.05% Thimerasol
Storage Conditions:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.
Storage Stability:	>2 years

Applications

Product Use:	Western blotting Antibody microarray
Antibody Dilution Recommended:	2 µg/ml for immunoblotting
Antibody Species Reactivity:	Human Chimpanzee Rhesus Macaque Dog

This product is for *in vitro* research use only and is not intended for use in humans or animals.

For more information on our products please visit <https://kinexus-ca.myshopify.com/> or contact us at 1-866-546-3987