

AB-NK085-4

JAK2-4 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-tyrosine kinase JAK2



KINEXUS

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Target Protein

Name Long:	Janus protein-tyrosine kinase 2
Alias:	JAK-2; Janus kinase 2; JTK10; Kinase Jak2 ; JAK2_ENST00000381652; JAK2_ENST00000381652; CCDS6457; A8K910; A9QVL8; B4DYV1; Q506Q0; Q52R86; Q8IXP2; ENSG00000096968
UniProt ID:	O60674
Sequence Predicted Mass (KDa):	130.674 (1132 AA; O60674)
Observed SDS-PAGE Mass (KDa):	125-135

Immunogen

Antibody Immunogen Source:	Human JAK2 sequence peptide Cat. No.: PE-01BBV95
Antibody Immunogen Sequence:	ELLTENDMLPNMRIGC
Location in Target:	Corresponds to amino acid residues E814 to G828; Post-kinase last alpha-chain
Peptide Type:	For pan-specific recognition of target expression levels.
Target Phosphosite:	Not phosphorylated

Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipепptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
Antibody Amount:	25 µg
Antibody Concentration:	1 mg/ml
Lot Number:	140724
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

AB-NK085-4

JAK2-4 Antibody



KINEXUS

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Applications

Product Use:	Western blotting Antibody microarrays
Antibody Dilution Recommended:	2 µg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Antibody Positive Controls:	Medium-strong immunoreactivity with recombinant human JAK2 on protein dot blots.
Overall Antibody Specificity:	High selectivity
Antibody Cross Reactivities:	No immunoreactivity on protein dot blots with recombinant human JAK1, JAK3 and Tyk2.

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-866-KINEXUS(546-3987)