## AB-NK092-2 LATS2-CT Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-serine/threonine kinase LATS2 (KPM)



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Target Protein	
Name Long:	Large tumour suppressor 2 protein-serine kinase
Alias:	FLJ13161; KPM; ENSG00000150457
UniProt ID:	Q9NRM7
Sequence Predicted Mass (KDa):	120.136 (1088 AA; Q9NRM7)
Observed SDS-PAGE Mass (KDa):	120-145
Immunogen	
	Human LATS2 (KRM) acquiance pontide Cat. No : DE 010005
Antibody Immunogen Source:	Human LATS2 (KPM) sequence peptide Cat. No.: PE-01BDO95 SQAESSDLESSDLVDQTEGC
Antibody Immunogen Sequence: Location in Target:	
	Corresponds to amino acid residues S1064 to C1083; C-terminus
Peptide Type:	For pan-specific recognition of target expression levels.
Target Phosphosite:	Not phosphorylated
Production	
Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the
	immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
Antibody Amount:	25 µg
Antibody Concentration:	0.5 mg/ml
Lot Number:	160807
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.





Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications	
Product Use:	Western blotting   Antibody microarrays
Antibody Dilution Recommended:	2 µg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals

This product is for in vitro research use only and is not intended for use in humans or animals.