AB-NK099-2 MEK1-CT Antibody

Pan-specific polyclonal antibody for monitoring the expression of human dual specificity protein kinase MEK1 (MKK1, MAP2K1)



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

get	

Name Long:	MAPK/ERK protein-serine kinase 1 (MKK1); Dual specificity mitogen-activated protein kinase kinase 1	
Alias:	ERK activator kinase 1; Kinase MEK1; MAP kinase kinase 1; MAP2K1; MAPK,ERK kinase 1; MAPK/ERK kinase 1; MAPKK 1; MAPKK1; MKK1; PRKMK1; CCDS10216.1; Q02750; ENSG00000169032	
UniProt ID:	Q02750	
Sequence Predicted Mass (KDa):	43.439 (393 AA; Q02750); 40.764 (367 AA; Q02750-2)	
Observed SDS-PAGE Mass (KDa):	44-47	

Immunogen

Antibody Immunogen Source:	Human MEK1 (MKK1, MAP2K1) sequence peptide Cat. No.: PE-01AQL90	
Antibody Immunogen Sequence:	HSFIKQSELEEVDFAGWLC	
Location in Target:	Corresponds to amino acid residues H358 to L375; C-terminus	
Peptide Type:	For pan-specific recognition of target expression levels.	
Target Phosphosite:	Not phosphorylated	

Production

Antibody Host Species:	Rabbit	
Antibody Type:	Polyclonal	
Antibody Ig Isotype Clone Lot:	Immunoglobulin G	
Production Method:	The immunizing peptide was produced on an A.B.I. automated peptide synthesizer. The peptide was cleaved from the synthesis resin with HF and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by amino acid analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.	
Antibody Amount:	25 μg	
Antibody Concentration:	0.5 mg/ml	
Lot Number:	10320	
Storage Buffer:	100 mM Tris-glycine, pH 7.0	
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months, but either 0.1% sodium azide or 0.05% Thimerasol should be added. Avoid repeated freeze-thaw cycles.	

AB-NK099-2 MEK1-CT Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications		
Product Use:	Western blotting Antibody microarrays	
Antibody Dilution Recommended:	2 μg/ml for immunoblotting	
Antibody Species Reactivity:	Human, mouse, rat and many other mammals as well as frog	
Antibody Positive Controls:	Strong immunoreactivity on Western blots with recombinant mouse MEK1 expressed in bacteria.	
Detection by Immunoblotting in Cell/Tissue Lysates:	With sea star oocyte extracts, it strongly immunoreacted with an appropriately sized 45 KDa protein that was also detected with STE7-VIII antibody, which al recognizes MEK1.	

This product is for in vitro research use only and is not intended for use in humans or animals.