AB-NK099-8 MEK1-1 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human dual specificity protein kinase MEK1 (MKK1, MAP2K1)



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

	1 -		
loro	Ot L	roti	מור
Targ			
	O	100	2111

Name Long:	MAPK/ERK protein-serine kinase 1 (MKK1); Dual specificity mitogen-activated	
	protein kinase kinase 1	
Alias:	ERK activator kinase 1; Kinase MEK1; MAP kinase kinase 1; MAP2K1; MAPK,ERK kinase 1; MAPK/ERK kinase 1; MAPKK 1; MAPKK1; MKK1; PRKMK1; CCDS10216.1; Q02750; ENSG00000169032	
UniProt ID:	Q02750	
Sequence Predicted Mass (KDa):	43.439 (393 AA; Q02750); 40.764 (367 AA; Q02750-2)	
Observed SDS-PAGE Mass (KDa):	: 44-47	

Immunogen

Antibody Immunogen Source:	Human MEK1 (MKK1, MAP2K1) sequence peptide Cat. No.: PE-01AVS80	
Antibody Immunogen Sequence:	PSGVFSLEFQDFVNKC	
Location in Target:	Corresponds to amino acid residues P326 to C341;	
Peptide Type:	For pan-specific recognition of target expression levels.	
Target Phosphosite:	Not phosphorylated	

Production

Antibody Host Species:	Rabbit	
Antibody Type:	Polyclonal	
Antibody Ig Isotype Clone Lot:	Immunoglobulin G	
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.	
Antibody Amount:	25 μg	
Antibody Concentration:	1 mg/ml	
Lot Number:	140724	
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol	
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.	

AB-NK099-8 MEK1-1 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

-//	\sim	\sim 1	-	ТΙС	ns
	v	\sim 1	-	\mathbf{u}	

7 ippiroditorio	
Product Use:	Western blotting Antibody microarrays
Antibody Dilution Recommended: 2 µg/ml for immunoblotting	
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Antibody Positive Controls:	Strong immunoreactivity with recombinant human MEK1 on protein dot blots.
Overall Antibody Specificity:	High selectivity
Antibody Cross Reactivities:	Strong immunoreactivity on protein dot blots with recombinant human MEK2, and no immunoreactivity with MKK4 and MKK6.

This product is for in vitro research use only and is not intended for use in humans or animals.