AB-NK101-6 MKK3-3 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human dual specificity protein kinase MKK3 (MAP2K3, MEK3)



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

	1		
101	$\boldsymbol{\alpha}$	_	מוח
101			ein
	\mathbf{g}		 . Ош

Name Long:	MAPK/ERK protein-serine kinase 3 beta isoform (MKK3 beta); Dual specificity
	mitogen-activated protein kinase kinase 3
Alias:	MAP kinase kinase 3; MAPK,ERK kinase 3; MAPK/ERK kinase 3; MAPKK 3; MAP2K3; MP2K3; PRKMK3; P46734; ENSG00000034152
UniProt ID:	P46734
Sequence Predicted Mass (KDa):	39.940 (352 AA; P46734-3); 39.318 (347 AA; P46734); 36.173 (318 AA; P46734-2)
Observed SDS-PAGE Mass (KDa):	37-40

Immunogen

Antibody Immunogen Source:	Human MKK3 (MAP2K3, MEK3) sequence peptide Cat. No.: PE-01AWD99
Antibody Immunogen Sequence:	CTDIAAFVKEILGEDS
Location in Target:	Corresponds to amino acid residues T333 to S347; C-terminus
Peptide Type:	For pan-specific recognition of target expression levels.
Target Phosphosite:	Not phosphorylated

Production

Antibody Host Species:	Rabbit	
Antibody Type:	Polyclonal	
Antibody Ig Isotype Clone Lot:	Immunoglobulin G	
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.	
Antibody Amount:	25 μg	
Antibody Concentration:	0.75 mg/ml	
Lot Number:	140724	
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol	
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.	

AB-NK101-6 MKK3-3 Antibody

Applications

Overall Antibody Specificity:

Antibody Cross Reactivities:



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

No immunoreactivity on protein dot blots with recombinant human MEK1, MKK4

Applications	
Product Use:	Western blotting Antibody microarrays
Antibody Dilution Recommended:	1 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals

Medium selectivity

and MKK6.

This product is for in vitro research use only and is not intended for use in humans or animals.