AB-NK104-5 MEK5-3 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human dual specificity protein kinase MEK5 (MAP2K5, MKK5)



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Target Protein

MAPK/ERK protein-serine kinase 5 (MKK5); Dual specificity mitogen-activated protein kinase kinase 5 (CXA1; HsT17454; MAP kinase kinase 5; MAP2K5; MAPK,ERK kinase 5; MAPK/ERK kinase 5; MAPKK5; MKK5; PRKMK5; HsT17454; ENSG00000137764

UniProt ID: Q13163

50.112 (448 AA; Q13163); 49.497 (444 AA; Q13163-3); 48.968 (438 AA;

Q13163-2); 46.131 (412 AA; Q13163-4)

Observed SDS-PAGE Mass (KDa): 47-52

Sequence Predicted Mass (KDa):

Immunogen

Antibody Immunogen Source: Human MEK5 (MAP2K5, MKK5) sequence peptide Cat. No.: PE-01AWA99

Antibody Immunogen Sequence: CRALEERRSQQGPP

Location in Target: Corresponds to amino acid residues C435 to P448; C-terminus

Peptide Type: For pan-specific recognition of target expression levels.

Target Phosphosite: Not phosphorylated

Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
Antibody Amount:	25 μg
Antibody Concentration:	0.5 mg/ml
Lot Number:	140724
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

AB-NK104-5 MEK5-3 Antibody

Applications

Antibody Cross Reactivities:



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

No immunoreactivity on protein dot blots with recombinant human MEK1, MKK4

Applications	
Product Use:	Western blotting Antibody microarrays
Antibody Dilution Recommended:	1 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals

and MKK6.

This product is for in vitro research use only and is not intended for use in humans or animals.