AB-NK108-4 MEKK2-2 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-serine/threonine kinase MEKK2 (MAP3K2)



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

66-74

Email: info@kinexus.ca Phone: 604-323-2547

 $\mathbf{r} \mathbf{a} \mathbf{a}$	_	מומי
rge		

Name Long:MAPK/ERK kinase kinase 2; Mitogen-activated protein kinase kinase 2Alias:Kinase MEKK2; M3K2; MAP/ERK kinase kinase 2; MAP3K2; MAPK/ERK kinase kinase 2; MAPKKK2; MEK kinase 2; MEKK 2; MEKK2B; ENSG00000169967UniProt ID:Q9Y2U5Sequence Predicted Mass (KDa):69.741 (619 AA; Q9Y2U5)

Immunogen

Observed SDS-PAGE Mass (KDa):

Antibody Immunogen Source:	Human MEKK2 (MAP3K2) sequence peptide Cat. No.: PE-01AXU90
Antibody Immunogen Sequence:	CLVINGSTQATNLEP
Location in Target:	Corresponds to amino acid residues L119 to P132
Peptide Type:	For pan-specific recognition of target expression levels.
Target Phosphosite:	Not phosphorylated

Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
Antibody Amount:	25 μg
Antibody Concentration:	0.03 / 0.5 mg/ml
Lot Number:	140801 / 160520
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

AB-NK108-4 MEKK2-2 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications							
		ne	tio	ica	m	n	Λ
Applications		115	UU	IUa	וטו	NU	

Product Use:	Western blotting Antibody microarrays
Antibody Dilution Recommended:	2 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Overall Antibody Specificity:	Very high selectivity
Antibody Cross Reactivities:	No immunoreactivity on protein dot blots with recombinant human MEKK6 and YSK4.

This product is for in vitro research use only and is not intended for use in humans or animals.