## AB-NK124-2 PAK6-APBD Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-serine/threonine kinase PAK6



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

### **Target Protein**

Name Long:

Protein-serine/threonine kinase PAK 6

PAK5; KIAA1264; P21 protein (Cdc42/Rac)-activated kinase 7; P21-activated kinase 7; PAK 7; PAK-5; PAK7; PAK-7; MGC26232; RP5-1119D9\_3; ENSG00000101349; KIAA1264KIAA1264; P21 protein (Cdc42/Rac)-activated kinase 7; P21-activated kinase 7; PAK 7; PAK-5; PAK7; PAK-7; MGC26232; RP5-1119D9\_3; ENSG00000101349; KIAA1264

UniProt ID:

Q9NQU5

Sequence Predicted Mass (KDa):

74.869 (681 AA; Q9NQU5); 69.777 (636 AA; Q9NQU5-2)

Observed SDS-PAGE Mass (KDa):

#### Immunogen

Antibody Immunogen Source: Human PAK6 sequence peptide Cat. No.: PE-01BEH90

Antibody Immunogen Sequence: CGSAMPVDGYISGLLNDIQK

Location in Target: Corresponds to amino acid residues G73 to K91; Post-PBD

Peptide Type: For pan-specific recognition of target expression levels.

Target Phosphosite: Not phosphorylated

#### **Production**

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
Antibody Amount:	25 μg
Antibody Concentration:	1 mg/ml
Lot Number:	160630
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

# AB-NK124-2 PAK6-APBD Antibody

Applications

**Antibody Cross Reactivities:** 



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

No significant cross-reactive proteins detected in diverse mouse tissues, except for a ~65 KDa protein in mouse testes, and weaker 50 and 37 KDa proteins in

mouse brain, skeletal muscle and thymus. No significant cross-reactivities detected in various human cancer cell lines, including A431, HEK-293, HeLa, HepG2, Jurkat and MCF7 cells, except for a strong 38 KDa proten in A431 cells.

Product Use:	Western blotting   Antibody microarrays
<b>Antibody Dilution Recommended:</b>	2 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Detection by Immunoblotting in Cell/Tissue Lysates:	Very strong immunoreactivity of a target-sized protein by Western blotting in various mouse brain.
Overall Antibody Specificity:	High selectivity

This product is for in vitro research use only and is not intended for use in humans or animals.