## AB-NK130-9 Akt2-3 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-serine/threonine kinase Akt2 (PKBb)



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

|  | ge |  |  |  |
|--|----|--|--|--|
|  |    |  |  |  |
|  |    |  |  |  |
|  |    |  |  |  |

| Name Long:                     | RAC-beta serine/threonine-protein kinase  |  |  |
|--------------------------------|---|--|--|
| Alias:                         | Akt2; PKB beta; RAC-beta serine,threonine protein kinase; RAC-PK-beta; PKBB; PKKBB; HIHGHH; PKBBETA; RAC-BETA; CCDS12552.1; ENSG00000105221 |  |  |
| UniProt ID:                    | P31751  |  |  |
| Sequence Predicted Mass (KDa): | 55.769 (481 AA; P31751); 51.083 (438 AA; P31751-2)  |  |  |
| Observed SDS-PAGE Mass (KDa):  | 55-60   |  |  |

## Immunogen

| Antibody Immunogen Source:   | Human Akt2 (PKBb) sequence peptide Cat. No.: PE-01AWR95         |
|------------------------------|---|
| Antibody Immunogen Sequence: | CDRYDSLGLLELDQRT  |
| Location in Target:          | Corresponds to amino acid residues D454 to T468; Pre-C-terminus |
| Peptide Type:                | For pan-specific recognition of target expression levels.       |
| Target Phosphosite:          | Not phosphorylated  |

## Production

| Antibody Host Species:            | Rabbit   |  |  |
|-----------------------------------|--|--|--|
| Antibody Type:                    | Polyclonal   |  |  |
| Antibody Ig Isotype Clone Lot:    | Immunoglobulin G   |  |  |
| Production Method:                | The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris. |  |  |
| Antibody Amount:                  | 25 μg  |  |  |
| Antibody Concentration:           | 1 mg/ml  |  |  |
| Lot Number:                       | 130906 / 140724  |  |  |
| Storage Buffer:                   | Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol  |  |  |
| Storage Conditions and Stability: | For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.  |  |  |

## AB-NK130-9 Akt2-3 Antibody

Applications

**Antibody Positive Controls:** 

**Overall Antibody Specificity:** 

**Antibody Cross Reactivities:** 



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Medium immunoreactivity with recombinant human Akt2 on protein dot blots.

No immunoreactivity on protein dot blots with recombinant human Akt1 and Akt3.

| Applications                          |  |
|---------------------------------------|--|
| Product Use:                          | Western blotting   Antibody microarrays  |
| <b>Antibody Dilution Recommended:</b> | 2 μg/ml for immunoblotting               |
| <b>Antibody Species Reactivity:</b>   | Human, mouse, rat and many other mammals |

High selectivity

This product is for in vitro research use only and is not intended for use in humans or animals.