AB-NK156-6 BRaf-2 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-serine/threonine kinase Raf-B (BRaf)



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Target Protein

Name Long:	B-Raf proto-oncogene serine-threonine-protein kinase; RafB proto-oncogene-
	encoded protein-serine kinase
Alias:	B-RAF; B-RAF proto-oncogene serine/threonine-protein kinase; BRAF1; C-RMIL; Kinase B-Raf; P94; RAFB1; RMIL; RMIL serine/threonine-protein kinase; V-Raf murine sarcoma viral oncogene homolog B1; MGC126806; MGC138284; CCDS5863.1; ENSG00000157764
UniProt ID:	P15056
Sequence Predicted Mass (KDa):	84.437 (766 AA; P15056)
Observed SDS-BAGE Mass (KDa):	87-07

Immunogen

Antibody Immunogen Source:	Human Raf-B (BRaf) sequence peptide Cat. No.: PE-01ASL95
Antibody Immunogen Sequence:	CGEKKPIGWDTDISWL
Location in Target:	Corresponds to amino acid residues G203 to L217
Peptide Type:	For pan-specific recognition of target expression levels.
Target Phosphosite:	Not phosphorylated

Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
Antibody Amount:	25 μg
Antibody Concentration:	1 mg/ml
Lot Number:	140724
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

AB-NK156-6 BRaf-2 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications

F F S S S S S	
Product Use:	Western blotting Antibody microarrays
Antibody Dilution Recommended:	2 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Overall Antibody Specificity:	Medium selectivity
Antibody Cross Reactivities:	No immunoreactivity on protein dot blots with recombinant human Raf1.

This product is for in vitro research use only and is not intended for use in humans or animals.