# AB-NK161-3

### Ron-2 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-tyrosine kinase Ron (RONa)



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Targ		

Name Long:	Macrophage-stimulating protein receptor alpha chain	
Alias:	CD136; CD136 antigen; CDW136; C-met-related tyrosine kinase; Kinase RON; Macrophage stimulating 1 receptor; Macrophage-stimulating protein receptor precursor; MSP receptor; MST1R; PTK8; STK; CCDS2807.1; ENSG00000164078	
UniProt ID:	Q04912	
Sequence Predicted Mass (KDa):	152.271 (1400 AA; Q04912); 147.277 (1351 AA; Q04912-2); 97.315 (907 AA; Q04912-5); 69.265 (647 AA; Q04912-6); 57.721 (541 AA; Q04912-4); 52.519 (495 AA; Q04912-3)	
Observed SDS-PAGE Mass (KDa):	150-220	

#### Immunogen

Antibody Immunogen Source:	Human Ron (RONa) sequence peptide Cat. No.: PE-01AVE90
Antibody Immunogen Sequence:	CAVRPTFRVLVGEVEQ
Location in Target:	Corresponds to amino acid residues A1329 to Q1343
Peptide Type:	For pan-specific recognition of target expression levels.
Target Phosphosite:	Not phosphorylated

### Production

Antibody Host Species:	Rabbit	
Antibody Type:	Polyclonal	
Antibody Ig Isotype Clone Lot:	Immunoglobulin G	
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.	
Antibody Amount:	25 μg	
Antibody Concentration:	1 mg/ml	
Lot Number:	140724	
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol	
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.	

## AB-NK161-3 Ron-2 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

#### Applications

Product Use:	Western blotting   Antibody microarrays	
Antibody Dilution Recommended:	2 μg/ml for immunoblotting	
Antibody Species Reactivity:	Human, mouse, rat and many other mammals	
Antibody Cross Reactivities:	No immunoreactivity on protein dot blots with recombinant human Met.	

This product is for in vitro research use only and is not intended for use in humans or animals.