## AB-NK163-4 Ros-3 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-tyrosine kinase Ros



Email: info@kinexus.ca Phone: 604-323-2547

Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Target Protein	
Name Long:	Orosomucoid 1 receptor-tyrosine kinase; Proto-oncogene tyrosine-protein kinase
Alias:	ROS1 c-Ros oncogene 1, receptor tyrosine kinase; c-Ros-1; KROS; MCF3; Proto- oncogene tyrosine-protein kinase ROS; ROS1; V-ros UR2 sarcoma virus oncogene homolog 1 (avian); V-ros UR2 sarcoma virus oncogene homologue 1 (avian); MGC163394; ENSG00000169071
UniProt ID:	P08922
Sequence Predicted Mass (KDa):	263.915 (2347 AA; P08922)
Observed SDS-PAGE Mass (KDa):	200-300
Immunogen	
Antibody Immunogen Source:	Human Ros sequence peptide Cat. No.: PE-01AVI90
Antibody Immunogen Sequence:	CRPTFHRIQDQLQLFR
Location in Target:	Corresponds to amino acid residues R2205 to R2219; Kinase last alpha-chain
Peptide Type:	For pan-specific recognition of target expression levels.
Target Phosphosite:	Not phosphorylated
Production	
Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
	The immunizing peptide was produced by solid phase synthesis on a multipep
Production Method:	peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
Production Method: Antibody Amount:	was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH
	was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
Antibody Amount:	<ul> <li>was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.</li> <li>25 µg</li> </ul>
Antibody Amount: Antibody Concentration:	<ul> <li>was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.</li> <li>25 μg</li> <li>1 mg/ml</li> </ul>





Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications	
Product Use:	Western blotting   Antibody microarrays
Antibody Dilution Recommended:	2 µg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Antibody Positive Controls:	Strong immunoreactivity with recombinant human Ros on protein dot blots.
<b>Overall Antibody Specificity:</b>	High selectivity
Antibody Cross Reactivities:	No immunoreactivity on protein dot blots with recombinant human Ret.

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit <u>www.kinexusproducts.ca</u> or contact us at 1-866-KINEXUS(546-3987)