

# AB-NK189-4

## JNK2-3 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-serine/threonine kinase JNK2 (MAPK9)



# KINEXUS

Address: 8755 Ash Street, Suite 1  
Vancouver, British Columbia,  
Canada V6P 6T3

Email: [info@kinexus.ca](mailto:info@kinexus.ca)  
Phone: 604-323-2547

### Target Protein

<b>Name Long:</b>	Jun N-terminus protein-serine kinase (stress-activated protein kinase (SAPK)) 2
<b>Alias:</b>	c-Jun N-terminal kinase 2; JNK2alpha; JNK2-alpha-2; JNK2A; JNK2B; JNK2beta; JNK-55; Jun kinase; Kinase JNK2; MAPK9; MK09; p54a; PRKM9; Stress-activated protein kinase JNK2; JNK-55; SAPK1a; JNK2BETA; p54aSAPK; CCDS4453.1; ENSG00000050748
<b>UniProt ID:</b>	P45984
<b>Sequence Predicted Mass (KDa):</b>	48.311 (424 AA; P45984-4); 48.139 (424 AA; P45984); 44.223 (382 AA; P45984-3); 44.051 (382 AA; P45984-2); 27.334 (242 AA; P45984-5)
<b>Observed SDS-PAGE Mass (KDa):</b>	48-54

### Immunogen

<b>Antibody Immunogen Source:</b>	Human JNK2 (MAPK9) sequence peptide Cat. No.: PE-01AXP90
<b>Antibody Immunogen Sequence:</b>	MSDSKGD SQFY SVQVC
<b>Location in Target:</b>	Corresponds to amino acid residues M1 to V15; N-terminus
<b>Peptide Type:</b>	For pan-specific recognition of target expression levels.
<b>Target Phosphosite:</b>	Not phosphorylated

### Production

<b>Antibody Host Species:</b>	Rabbit
<b>Antibody Type:</b>	Polyclonal
<b>Antibody Ig Isotype Clone Lot:</b>	Immunoglobulin G
<b>Production Method:</b>	The immunizing peptide was produced by solid phase synthesis on a multipепptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
<b>Antibody Amount:</b>	25 µg
<b>Antibody Concentration:</b>	0.24 / 0.25 mg/ml
<b>Lot Number:</b>	140801 / 160520
<b>Storage Buffer:</b>	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
<b>Storage Conditions and Stability:</b>	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

# AB-NK189-4

## JNK2-3 Antibody



# KINEXUS

Address: 8755 Ash Street, Suite 1  
Vancouver, British Columbia,  
Canada V6P 6T3

Email: [info@kinexus.ca](mailto:info@kinexus.ca)  
Phone: 604-323-2547

### Applications

<b>Product Use:</b>	Western blotting   Antibody microarrays
<b>Antibody Dilution Recommended:</b>	2 µg/ml for immunoblotting
<b>Antibody Species Reactivity:</b>	Human, mouse, rat and many other mammals
<b>Overall Antibody Specificity:</b>	Very high selectivity
<b>Antibody Cross Reactivities:</b>	No immunoreactivity on protein dot blots with recombinant human JNK1 and JNK3.

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit [www.kinexusproducts.ca](http://www.kinexusproducts.ca) or contact us at 1-866-KINEXUS(546-3987)