## AB-NK197-2 JNK3-1 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-serine/threonine kinase JNK3 (MAPK10)



Email: info@kinexus.ca Phone: 604-323-2547

Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Target Protein	
Name Long:	Jun N-terminus protein-serine kinase (stress-activated protein kinase (SAPKb)) 3
Alias:	c-Jun N-terminal kinase 3; FLJ12099; FLJ33785; JNK3A; JNK3-alpha-2; Kinase JNK3; MAP kinase p49 3F12; MAPK10; MGC50974; MK10; P493F12; p54bSAPK; PRKM10; Stress-activated protein kinase JNK3; PRKM10; CCDS3613.1; ENSG00000109339
UniProt ID:	P53779
Sequence Predicted Mass (KDa):	52.585 (464 AA; P53779); 48.554 (422 AA; P53779-2); 48.128 (426 AA; P53779-3); 31.933 (277 AA; P53779-4)
Observed SDS-PAGE Mass (KDa):	48-54
Immunogen	
Antibody Immunogen Source:	Human JNK3 (MAPK10) sequence peptide Cat. No.: PE-01AXQ99
Antibody Immunogen Sequence:	CAGLTFPKLFPDSLFP
Location in Target:	Corresponds to amino acid residues A305 to P319; Pre-kinase last alpha-chain
Peptide Type:	For pan-specific recognition of target expression levels.
Target Phosphosite:	Not phosphorylated
Production Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Type. Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
Antibody Amount:	25 μg
Antibody Concentration:	1 mg/ml
Lot Number:	140724
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.





Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications	
Product Use:	Western blotting   Antibody microarrays
Antibody Dilution Recommended:	2 µg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
<b>Overall Antibody Specificity:</b>	Very high selectivity
Antibody Cross Reactivities:	No immunoreactivity on protein dot blots with recombinant human JNK1 and JNK2.

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit <u>www.kinexusproducts.ca</u> or contact us at 1-866-KINEXUS(546-3987)