AB-NK225-3P MEKK6 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-serine/threonine kinase MEKK6 (MAP3K6)



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

| Target Protein | |
|---|--|
| Name Long: | Mitogen-activated protein kinase kinase kinase 6 |
| Alias: | Apoptosis signal regulating kinase 2; Kinase MAP3K6; M3K6; MAP3K6; MAPKKK6; MEKK6; MEK kinase 6; MGC125653; MGC20114; MGC125653; CCDS299.1; ENSG00000142733 |
| UniProt ID: | O95382 |
| Sequence Predicted Mass (KDa): | 142.596 (1288 AA; O95382); 141.648 (1280 AA; O95382-3); 88.674 (797 AA; O95382-2) |
| Observed SDS-PAGE Mass (KDa): | 110-130 |
| Immunogen | |
| Antibody Immunogen Source: | Human MEKK6 (MAP3K6) sequence peptide Cat. No.: PE-01ASS95 |
| Antibody Immunogen Sequence: | CRIWRAILAQRAGSTP |
| Location in Target: | Corresponds to amino acid residues C1268 to P1283; Pre-C-terminus |
| Peptide Type: | For pan-specific recognition of target expression levels. |
| Target Phosphosite: | Not phosphorylated |
| Production | |
| Antibody Host Species: | Rabbit |
| | Rabbit Polyclonal |
| Antibody Host Species: | |
| Antibody Host Species: Antibody Type: | Polyclonal |
| Antibody Host Species: Antibody Type: Antibody Ig Isotype Clone Lot: | Polyclonal Immunoglobulin G The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH |
| Antibody Host Species: Antibody Type: Antibody Ig Isotype Clone Lot: Production Method: | Polyclonal Immunoglobulin G The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris. |
| Antibody Host Species: Antibody Type: Antibody Ig Isotype Clone Lot: Production Method: Antibody Amount: | Polyclonal Immunoglobulin G The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris. 25 μg |
| Antibody Host Species: Antibody Type: Antibody Ig Isotype Clone Lot: Production Method: Antibody Amount: Antibody Concentration: | PolyclonalImmunoglobulin GThe immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.25 μg1 mg/ml |

AB-NK225-3P MEKK6 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

| Applications | |
|--------------------------------|---|
| Product Use: | Western blotting Antibody microarrays |
| Antibody Dilution Recommended: | 2 µg/ml for immunoblotting |
| Antibody Species Reactivity: | Human, mouse, rat and many other mammals |
| Overall Antibody Specificity: | Very high selectivity |
| Antibody Cross Reactivities: | No immunoreactivity on protein dot blots with recombinant human MEKK2 and YSK4. |

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit <u>www.kinexusproducts.ca</u> or contact us at 1-866-KINEXUS(546-3987)