

AB-NK226-2P

VEGFR1 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-tyrosine kinase VEGFR1 (Flt1)



KINEXUS

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Target Protein

| | |
|---------------------------------------|--|
| Name Long: | Vascular endothelial growth factor receptor 1 |
| Alias: | FLT; Flt-1; Fms-like tyrosine kinase 1; Fms-related tyrosine kinase 1; FRT; Tyrosine-protein kinase FRT; Tyrosine-protein kinase receptor FLT; VGFR1; CCDS9330.1; ENSG00000102755 |
| UniProt ID: | P17948 |
| Sequence Predicted Mass (KDa): | 150.769 (1338 AA; P17948); 82.124 (733 AA; P17948-3); 77.474 (687 AA; P17948-2); 62.954 (556 AA; P17948-5); 60.917 (541 AA; P17948-4); 52.613 (463 AA; P17948-6); 41.175 (361 AA; P17948-8); 39.148 (343 AA; P17948-7) |
| Observed SDS-PAGE Mass (KDa): | 150-190 |

Immunogen

| | |
|-------------------------------------|--|
| Antibody Immunogen Source: | Human VEGFR1 (Flt1) sequence peptide Cat. No.: PE-01ASD85 |
| Antibody Immunogen Sequence: | PEIIHMTEGRELVIPIC |
| Location in Target: | Corresponds to amino acid residues P143 to C158; Just before and part of the second Ig-like C2-type domain |
| Peptide Type: | For pan-specific recognition of target expression levels. |
| Target Phosphosite: | Not phosphorylated |

Production

| | |
|--|---|
| Antibody Host Species: | Rabbit |
| Antibody Type: | Polyclonal |
| Antibody Ig Isotype Clone Lot: | Immunoglobulin G |
| Production Method: | The immunizing peptide was produced by solid phase synthesis on a multipепptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris. |
| Antibody Amount: | 25 µg |
| Antibody Concentration: | mg/ml |
| Lot Number: | |
| Storage Buffer: | Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol |
| Storage Conditions and Stability: | For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles. |

AB-NK226-2P

VEGFR1 Antibody



KINEXUS

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Applications

| | |
|---------------------------------------|--|
| Product Use: | Western blotting Antibody microarrays |
| Antibody Dilution Recommended: | 2 µg/ml for immunoblotting |
| Antibody Species Reactivity: | Human, mouse, rat and many other mammals |
| Antibody Positive Controls: | Very strong immunoreactivity with recombinant human VEGFR1 on protein dot blots. |
| Overall Antibody Specificity: | High selectivity |
| Antibody Cross Reactivities: | No immunoreactivity on protein dot blots with recombinant human VGFR2 and VGFR3. |

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-866-KINEXUS(546-3987)