## AB-NK229-3 **MELK-2** Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-serine/threonine kinase MELK



Email: info@kinexus.ca

Phone: 604-323-2547

Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Target Protein	
Name Long:	Maternal embryonic leucine zipper kinase
Alias:	hMELK; hPK38; KIAA0175; Maternal embryonic leucine zipper kinase; MELK; pEg3 kinase; Protein kinase PK38
UniProt ID:	Q14680
Sequence Predicted Mass (KDa):	74.642 (651 AA; Q14680); 71.174 (619 AA; Q14680-6); 70.150 (610 AA; Q14680-7); 69.116 (603 AA; Q14680-8); 66.547 (580 AA; Q14680-5); 66.399 (580 AA; Q14680-2); 59.576 (520 AA; Q14680-4); 52.528 (457 AA; Q14680-3)
Observed SDS-PAGE Mass (KDa):	70-75
Immunogen	
Antibody Immunogen Source:	Human MELK sequence peptide Cat. No.: PE-01ASU90

Antibody Immunogen Source:	Human MELK sequence peptide Cat. No.: PE-01ASU90
Antibody Immunogen Sequence:	CKKRISMKNLLNHPWI
Location in Target:	Corresponds to amino acid residues K249 to I263; Kinase last alpha-chain
Peptide Type:	For pan-specific recognition of target expression levels.
Target Phosphosite:	Not phosphorylated

Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
Antibody Amount:	25 µg
Lot Number:	140724
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.





Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications	
Product Use:	Western blotting   Antibody microarrays
Antibody Dilution Recommended:	2 µg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Antibody Positive Controls:	Medium-strong immunoreactivity with recombinant human MELK on protein dot blots.
<b>Overall Antibody Specificity:</b>	High selectivity

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit <u>www.kinexusproducts.ca</u> or contact us at 1-866-KINEXUS(546-3987)