AB-NK231-2P ErbB3 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-tyrosine kinase ErbB3 (HER3)



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Target Protein

| Name Long: | ErbB3 (HER3) receptor-tyrosine kinase |
|--------------------------------|--|
| Alias: | C-erbB3; ERB3; HER3; Kinase ErbB3; LCCS2; Receptor protein-tyrosine kinase erbB-3 precursor; Tyrosine kinase-type cell surface receptor HER3; V-erb-b2 erythroblastic leukemia viral oncogene 3; ErbB-3; c-erbB3; erbB3-S; MDA-BF-1; MGC88033; c-erbB-3; p180-ErbB3; p45-sErbB3; p85-sErbB3; ENSG00000065361; O75810; O75811; O75812; O75813; Q9BUD7; Q9NNX2; Q9NNX3 |
| UniProt ID: | P21860 |
| Sequence Predicted Mass (KDa): | 148.098 (1342 AA; P21860); 141.752 (1283 AA; P21860-4); 77.426 (699 AA; P21860-5); 36.490 (331 AA; P21860-3); 20.136 (183 AA; P21860-2) |
| Observed SDS-PAGE Mass (KDa): | 160-200 |

Immunogen

| Antibody Immunogen Source: | Human ErbB3 (HER3) sequence peptide Cat. No.: PE-01ARS99 |
|------------------------------|---|
| Antibody Immunogen Sequence: | DRDAEIVVKDNGRSC |
| Location in Target: | Corresponds to amino acid residues D169 to C183; |
| Peptide Type: | For pan-specific recognition of target expression levels. |
| Target Phosphosite: | Not phosphorylated |

Production

| Antibody Host Species: | Rabbit |
|-----------------------------------|--|
| Antibody Type: | Polyclonal |
| Antibody Ig Isotype Clone Lot: | Immunoglobulin G |
| Production Method: | The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris. |
| Antibody Amount: | 25 μg |
| Storage Buffer: | Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol |
| Storage Conditions and Stability: | For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles. |

AB-NK231-2P ErbB3 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications

| Product Use: | Western blotting Antibody microarrays |
|---------------------------------------|--|
| Antibody Dilution Recommended: | 2 μg/ml for immunoblotting |
| Antibody Species Reactivity: | Human, mouse, rat and many other mammals |
| Antibody Positive Controls: | Very strong immunoreactivity with recombinant human ErbB3 on protein dot blots. |
| Overall Antibody Specificity: | Medium selectivity |
| Antibody Cross Reactivities: | No immunoreactivity on protein dot blots with recombinant human ErbB2 and ErbB4. |

This product is for in vitro research use only and is not intended for use in humans or animals.