# AB-NK238-1 Abl2-1 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-tyrosine kinase Abl2 (Arg)



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

### **Target Protein**

| Name Long:                     | Abelson murine leukemia viral oncogene homologue 2   |
|--------------------------------|--|
| Alias:                         | Abelson-related gene; ABLL; ARG; Kinase Arg; Tyrosine kinase ARG; v-Abl Abelson murine leukemia viral oncogene 2; v-Abl Abelson murine leukemia viral oncogene homologue 2; hCG_23565; FLJ222224; FLJ31718; FLJ41441; RP11-177A2_3; ENSG00000143322; D1MPS6      |
| UniProt ID:                    | P42684   |
| Sequence Predicted Mass (KDa): | 128.343 (1182 AA; P42684); 126.684 (1167 AA; P42684-3); 126.219 (1161 AA; P42684-6); 124.560 (1146 AA; P42684-2); 117.942 (1079 AA; P42684-5); 116.283 (1064 AA; P42684-10); 115.817 (1058 AA; P42684-7); 114.159 (1043 AA; P42684-4); 60.864 (542 AA; P42684-8) |
| Observed SDS-PAGE Mass (KDa):  | 130-140  |

#### Immunogen

| Antibody Immunogen Source:   | Human Abl2 (Arg) sequence peptide Cat. No.: PE-01ATM99                   |
|------------------------------|--|
| Antibody Immunogen Sequence: | CDSISEEVAEELGRAA   |
| Location in Target:          | Corresponds to amino acid residues S547 to A560; Kinase last alpha-chain |
| Peptide Type:                | For pan-specific recognition of target expression levels.                |
| Target Phosphosite:          | Not phosphorylated   |

## Production

| Antibody Host Species:            | Rabbit   |
|-----------------------------------|--|
| Antibody Type:                    | Polyclonal   |
| Antibody Ig Isotype Clone Lot:    | Immunoglobulin G   |
| Production Method:                | The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris. |
| Antibody Amount:                  | 25 μg  |
| Antibody Concentration:           | 0.5 mg/ml  |
| Lot Number:                       | 140724   |
| Storage Buffer:                   | Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol  |
| Storage Conditions and Stability: | For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.  |

# AB-NK238-1 Abl2-1 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

### Applications

| F F S S S S S                         |   |
|---------------------------------------|---|
| Product Use:                          | Western blotting   Antibody microarrays                               |
| <b>Antibody Dilution Recommended:</b> | 2 μg/ml for immunoblotting  |
| Antibody Species Reactivity:          | Human, mouse, rat and many other mammals                              |
| <b>Antibody Cross Reactivities:</b>   | No immunoreactivity on protein dot blots with recombinant human Abl1. |

This product is for in vitro research use only and is not intended for use in humans or animals.