

# AB-NK248-1

## p38b-1 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-serine/threonine kinase p38b MAPK (MAPK11)



**KINEXUS**

Address: 8755 Ash Street, Suite 1  
Vancouver, British Columbia,  
Canada V6P 6T3

Email: [info@kinexus.ca](mailto:info@kinexus.ca)  
Phone: 604-323-2547

### Target Protein

<b>Name Long:</b>	Mitogen-activated protein-serine kinase p38 beta; Mitogen-activated protein kinase 11
<b>Alias:</b>	MAPK11; MK11; P38 MAPK-beta; P38-2; P38b; PRKM11; SAPK2; Stress-activated protein kinase-2; p28b MAPK; p38Beta; CCDS14090.1; Q15759; ENSG00000185386
<b>UniProt ID:</b>	Q15759
<b>Sequence Predicted Mass (KDa):</b>	41.357 (364 AA; Q15759); 23.603 (213 AA; Q15759-3)
<b>Observed SDS-PAGE Mass (KDa):</b>	40-48

### Immunogen

<b>Antibody Immunogen Source:</b>	Human p38b MAPK (MAPK11) sequence peptide Cat. No.: PE-01AYD85
<b>Antibody Immunogen Sequence:</b>	MSGPRAGFYRQELNKC
<b>Location in Target:</b>	Corresponds to amino acid residues M1 to K15; N-terminus
<b>Peptide Type:</b>	For pan-specific recognition of target expression levels.
<b>Target Phosphosite:</b>	Not phosphorylated

### Production

<b>Antibody Host Species:</b>	Rabbit
<b>Antibody Type:</b>	Polyclonal
<b>Antibody Ig Isotype Clone Lot:</b>	Immunoglobulin G
<b>Production Method:</b>	The immunizing peptide was produced by solid phase synthesis on a multipеп peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
<b>Antibody Amount:</b>	25 µg
<b>Antibody Concentration:</b>	1 mg/ml
<b>Lot Number:</b>	140724
<b>Storage Buffer:</b>	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
<b>Storage Conditions and Stability:</b>	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

# AB-NK248-1

## p38b-1 Antibody



# KINEXUS

Address: 8755 Ash Street, Suite 1  
Vancouver, British Columbia,  
Canada V6P 6T3

Email: [info@kinexus.ca](mailto:info@kinexus.ca)  
Phone: 604-323-2547

### Applications

<b>Product Use:</b>	Western blotting   Antibody microarrays
<b>Antibody Dilution Recommended:</b>	1 µg/ml for immunoblotting
<b>Antibody Species Reactivity:</b>	Human, mouse, rat and many other mammals
<b>Antibody Positive Controls:</b>	Strong immunoreactivity with recombinant human p38b on protein dot blots.
<b>Overall Antibody Specificity:</b>	High selectivity
<b>Antibody Cross Reactivities:</b>	Medium immunoreactivity on protein dot blots with recombinant human p38a and p38g, and no immunoreactivity on protein dot blots with recombinant human p38d.

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit [www.kinexusproducts.ca](http://www.kinexusproducts.ca) or contact us at 1-866-KINEXUS(546-3987)