# AB-NK253-1 WNK2-1 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-serine/threonine kinase WNK2



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

### **Target Protein**

Name Long:

Serine/threonine-protein kinase WNK2

Antigen NY-CO-43; KIAA1760; NY-CO-43; P,OKcl.13; P/OKcl.13; PRKWNK2;

Protein kinase with no lysine 2; Protein kinase, lysine-deficient 2; WNK lysine deficient protein kinase 2; SDCCAG43; KAIA302; ENSG00000165238

UniProt ID:

Q9Y3S1

242.676 (2297 AA; Q9Y3S1); 238.279 (2254 AA; Q9Y3S1-4); 234.948 (2224 AA;

Sequence Predicted Mass (KDa): 242.076 (2297 AA, Q913S1), 236.279 (2234 AA, Q913S1-2); 84.868 (779 AA; Q9Y3S1-3)

Observed SDS-PAGE Mass (KDa): 250-300

#### Immunogen

Antibody Immunogen Source: Human WNK2 sequence peptide Cat. No.: PE-01AZB99

Antibody Immunogen Sequence: MDGDGGRRDVPGTLMC

Location in Target: Corresponds to amino acid residues M1 to M15; N-terminus

Peptide Type: For pan-specific recognition of target expression levels.

Target Phosphosite: Not phosphorylated

#### Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
Antibody Amount:	25 μg
Antibody Concentration:	1 mg/ml
Lot Number:	140724
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

# AB-NK253-1 WNK2-1 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

### Applications

, ipprisations	
Product Use:	Western blotting   Antibody microarrays
<b>Antibody Dilution Recommended:</b>	2 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Overall Antibody Specificity:	High selectivity
<b>Antibody Cross Reactivities:</b>	No immunoreactivity on protein dot blots with recombinant human WNK1.

This product is for in vitro research use only and is not intended for use in humans or animals.