AB-NK256-1 YSK4-1 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-serine/threonine kinase YSK4 (MAP3K19)



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Target Protein	
Name Long:	SPS1/STE20-related protein kinase YSK4; Mitogen-activated protein kinase kinase kinase 8; MAP3K19
Alias:	FLJ23074; Rck; Regulated in COPD, protein kinase; Ysk4; YSK4 Sps1/Ste20- related kinase
UniProt ID:	Q56UN5
Sequence Predicted Mass (KDa):	150.537 (1328 AA; Q56UN5); 137.548 (1215 AA; Q56UN5-3); 127.941 (1127 AA; Q56UN5-2); 57.659 (510 AA; Q56UN5-5); 52.504 (462 AA; Q56UN5-4); 51.920 (460 AA; Q56UN5-7); 18.702 (168 AA; Q56UN5-6)
Observed SDS-PAGE Mass (KDa):	50-70
Immunogen	
Antibody Immunogen Source:	Human YSK4 (MAP3K19) sequence peptide Cat. No.: PE-01AZK85
Antibody Immunogen Sequence:	MSSMPKPERHAESLLC
Location in Target:	Corresponds to amino acid residues M1 to L15; N-terminus
Peptide Type:	For pan-specific recognition of target expression levels.
Target Phosphosite:	Not phosphorylated
Production	
Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Deciding Mathedu	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purit was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously

Production Method:

immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.

Antibody Amount:	25 μg
Antibody Concentration:	0.5 mg/ml
Lot Number:	140724
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.





Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications	
Product Use:	Western blotting Antibody microarrays
Antibody Dilution Recommended:	2 µg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Overall Antibody Specificity:	High selectivity
Antibody Cross Reactivities:	No immunoreactivity on protein dot blots with recombinant human MEKK2 and MEKK6.

This product is for in vitro research use only and is not intended for use in humans or animals.