## AB-NK257-2 BCKD-2 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-serine/threonine kinase BCKD (BCKDK)



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

#### **Target Protein**

Name Long: [3-methyl-2-oxobutanoate dehydrogenase [lipoamide]] kinase, mitochondrial

Alias: BCKD; BCKDHKIN; BCKD-kinase; Branched-chain alpha-ketoacid

dehydrogenase kinase; EC 2.7.1.115; EC 2.7.11.4

UniProt ID: O14874

Sequence Predicted Mass (KDa): 46.360 (412 AA; O14874); 41.090 (365 AA; O14874-3); 37.768 (335 AA;

O14874-2)

Observed SDS-PAGE Mass (KDa): 45-50

#### Immunogen

Antibody Immunogen Source: Human BCKD (BCKDK) sequence peptide Cat. No.: PE-01AZU90

Antibody Immunogen Sequence: CILHVHELYIRAFQKL

Location in Target: Corresponds to amino acid residues I122 to L136;

**Peptide Type:** For pan-specific recognition of target expression levels.

Target Phosphosite: Not phosphorylated

#### Production

Antibody Host Species: Rabbit
Antibody Type: Polyclonal

Antibody Ig Isotype Clone Lot: Immunoglobulin G

The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by

was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously

Production Method:

injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with

0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH

7.0 with saturated Tris.

Antibody Amount: 25 μg

Antibody Concentration: 1 mg/ml Lot Number: 140724

**Storage Buffer:** Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol

Storage Conditions and Stability:

For long term storage, keep frozen at -40°C or lower. Stock solution can be kept

at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

# AB-NK257-2 BCKD-2 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

### Applications

F P - F F F F F	
Product Use:	Western blotting   Antibody microarrays
<b>Antibody Dilution Recommended:</b>	2 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Overall Antibody Specificity:	High selectivity

This product is for in vitro research use only and is not intended for use in humans or animals.