## AB-NK265-1 **DMPK2-NT Antibody**

Pan-specific polyclonal antibody for monitoring the expression of human protein-serine/threonine kinase DMPK2 (MRCK-q)



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

#### **Target Protein**

Name Long: Myotonic dystrophy protein kinase 2; Serine-threonine-protein kinase MRCK

MRCKG; Myotonic dystrophy protein kinase like protein; CDC42BPG; Alias:

**HSMDPKIN UniProt ID:** Q6DT37

Sequence Predicted Mass (KDa): 172.459 (1551 AA; Q6DT37)

Observed SDS-PAGE Mass (KDa): n/a

#### **Immunogen**

Antibody Immunogen Source: Human DMPK2 (MRCK-g) sequence peptide Cat. No.: PE-01BCS80

**Antibody Immunogen Sequence:** MERRLRALEQLARGEAGGC

**Location in Target:** Corresponds to amino acid residues M1 to C19; N-terminus

Peptide Type: For pan-specific recognition of target expression levels.

**Target Phosphosite:** Not phosphorylated

#### Production

**Antibody Host Species:** Rabbit **Antibody Type:** Polyclonal

**Antibody Ig Isotype Clone Lot:** Immunoglobulin G

> The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity

was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously

**Production Method:** injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with

> 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.

25 µg **Antibody Amount:** 

**Antibody Concentration:** 0.5 mg/ml Lot Number: 160804

Storage Buffer: Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol

For long term storage, keep frozen at -40°C or lower. Stock solution can be kept Storage Conditions and Stability:

at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

# AB-NK265-1 DMPK2-NT Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

### Applications

· Pp.	
Product Use:	Western blotting   Antibody microarrays
Antibody Dilution Recommended:	2 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals

This product is for in vitro research use only and is not intended for use in humans or animals.