## AB-NK270-3 GSK3b-ANT Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-serine/threonine kinase GSK3b



Email: info@kinexus.ca

Phone: 604-323-2547

Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Target Protein	
Name Long:	Glycogen synthase-serine kinase 3 beta
Alias:	Factor A; Glycogen synthase kinase 3 beta; GSK-3 beta; GSK3-beta; Kinase GSK-3 beta; Kinase GSK3-beta; Protein kinase GSK-3-beta; CCDS2996.1; ENSG00000082701
UniProt ID:	P49841
Sequence Predicted Mass (KDa):	46.744 (420 AA; P49841-1); 48.034 (433 AA; P49841-2)
Observed SDS-PAGE Mass (KDa):	45-50
Immunogen	
Antibody Immunogen Source:	Human GSK3b sequence peptide Cat. No.: PE-01BDB85
Antibody Immunogen Sequence:	CKPVQQPSAFGSMKVSRDKD
Location in Target:	Corresponds to amino acid residues C14 to D33; Near N-terminus
Peptide Type:	For pan-specific recognition of target expression levels.
Target Phosphosite:	Not phosphorylated
Production	
Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
Antibody Amount:	25 μg
Antibody Concentration:	1 mg/ml
Lot Number:	160804
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

## AB-NK270-3 GSK3b-ANT Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications		
Product Use:	Western blotting   Antibody microarrays	
Antibody Dilution Recommended:	2 µg/ml for immunoblotting	
Antibody Species Reactivity:	Human, mouse, rat and many other mammals	

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit <u>www.kinexusproducts.ca</u> or contact us at 1-866-KINEXUS(546-3987)