

# AB-NK275-1

## MARK2-ANT Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-serine/threonine kinase MARK2



**KINEXUS**

Address: 8755 Ash Street, Suite 1  
Vancouver, British Columbia,  
Canada V6P 6T3

Email: [info@kinexus.ca](mailto:info@kinexus.ca)  
Phone: 604-323-2547

### Target Protein

<b>Name Long:</b>	MAP/microtubule affinity-regulating protein-serine kinase 2
<b>Alias:</b>	ELKL motif kinase 1; EMK; EMK1; Kinase MARK2; MAP/microtubule affinity-regulating kinase 2; MGC99619; PAR-1; Par1b; PAR-1b
<b>UniProt ID:</b>	Q7KZI7
<b>Sequence Predicted Mass (kDa):</b>	87.911 (788 AA; Q7KZI7); 86.762 (779 AA; Q7KZI7-8); 86.129 (773 AA; Q7KZI7-11); 84.333 (755 AA; Q7KZI7-6); 83.184 (746 AA; Q7KZI7-3); 83.113 (745 AA; Q7KZI7-14); 82.551 (740 AA; Q7KZI7-12); 82.429 (734 AA; Q7KZI7-9); 81.280 (725 AA; Q7KZI7-4); 81.209 (724 AA; Q7KZI7-16); 80.647 (719 AA; Q7KZI7-5); 79.427 (709 AA; Q7KZI7-15); 78.851 (701 AA; Q7KZI7-2); 77.702 (692 AA; Q7KZI7-7); 77.631 (691 AA; Q7KZI7-13); 77.069 (686 AA; Q7KZI7-10)
<b>Observed SDS-PAGE Mass (kDa):</b>	80-90

### Immunogen

<b>Antibody Immunogen Source:</b>	Human MARK2 sequence peptide Cat. No.: PE-01BDR80
<b>Antibody Immunogen Sequence:</b>	PTLGHLDSKPPSSKSNMIRGC
<b>Location in Target:</b>	Corresponds to amino acid residues P19 to G37; Post-N-terminus
<b>Peptide Type:</b>	For pan-specific recognition of target expression levels.
<b>Target Phosphosite:</b>	Not phosphorylated

### Production

<b>Antibody Host Species:</b>	Rabbit
<b>Antibody Type:</b>	Polyclonal
<b>Antibody Ig Isotype Clone Lot:</b>	Immunoglobulin G
<b>Production Method:</b>	The immunizing peptide was produced by solid phase synthesis on a multipепptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
<b>Antibody Amount:</b>	25 µg
<b>Antibody Concentration:</b>	0.5 mg/ml
<b>Lot Number:</b>	160807
<b>Storage Buffer:</b>	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
<b>Storage Conditions and Stability:</b>	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

# AB-NK275-1

## MARK2-ANT Antibody



# KINEXUS

Address: 8755 Ash Street, Suite 1  
Vancouver, British Columbia,  
Canada V6P 6T3

Email: [info@kinexus.ca](mailto:info@kinexus.ca)  
Phone: 604-323-2547

### Applications

<b>Product Use:</b>	Western blotting   Antibody microarrays
<b>Antibody Dilution Recommended:</b>	2 µg/ml for immunoblotting
<b>Antibody Species Reactivity:</b>	Human, mouse, rat and many other mammals

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit [www.kinexusproducts.ca](http://www.kinexusproducts.ca) or contact us at 1-866-KINEXUS(546-3987)