AB-NK276-2 MARK3-BCT Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-serine/threonine kinase MARK3



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Target Protein

Name Long:	MAP/microtubule affinity-regulating protein-serine kinase 3
Alias:	Cdc25C- associated protein kinase 1; Cdc25C-associated protein kinase 1; CTAK1; C-TAK1; ELKL motif kinase 2; KP78; PAR1A; Protein kinase STK10; Ser/Thr protein kinase PAR-1; Serine/threonine protein kinase p78; Serine/threonine-protein kinase p78
UniProt ID:	P27448
Sequence Predicted Mass (KDa):	85.441 (760 AA; P27448-7); 84.489 (753 AA; P27448); 84.015 (752 AA; P27448-2); 83.308 (744 AA; P27448-4); 81.499 (729 AA; P27448-3); 79.936 (713 AA; P27448-6); 74.151 (659 AA; P27448-8)
Observed SDS-PAGE Mass (KDa):	80-85

Immunogen

Antibody Immunogen Source:	Human MARK3 sequence peptide Cat. No.: PE-01BDU90
Antibody Immunogen Sequence:	RFIKRLPTEYERNGRYEGC
Location in Target:	Corresponds to amino acid residues R620 to G637; Pre-C-terminus
Peptide Type:	For pan-specific recognition of target expression levels.
Target Phosphosite:	Not phosphorylated

Production

Antibody Type: Antibody Ig Isotype Clone Lot: Immunoglobulin G The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris. Antibody Amount: 25 µg Antibody Concentration: Lot Number: 160807 Storage Buffer: Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.	Antibody Host Species:	Rabbit
The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris. Antibody Amount: 25 μg Antibody Concentration: 0.5 mg/ml Lot Number: 160807 Storage Conditions and Stability: For long term storage, keep frozen at -40°C or lower. Stock solution can be kept	Antibody Type:	Polyclonal
peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris. Antibody Amount: 25 µg Antibody Concentration: 0.5 mg/ml Lot Number: 160807 Storage Buffer: Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol For long term storage, keep frozen at -40°C or lower. Stock solution can be kept	Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Antibody Concentration: Lot Number: 160807 Storage Buffer: Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol For long term storage, keep frozen at -40°C or lower. Stock solution can be kept	Production Method:	peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH
Lot Number: 160807 Storage Buffer: Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol Storage Conditions and Stability: For long term storage, keep frozen at -40°C or lower. Stock solution can be kept	Antibody Amount:	25 μg
Storage Buffer: Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol Storage Conditions and Stability: For long term storage, keep frozen at -40°C or lower. Stock solution can be kept	Antibody Concentration:	0.5 mg/ml
Storage Conditions and Stability: For long term storage, keep frozen at -40°C or lower. Stock solution can be kept	Lot Number:	160807
	Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
	Storage Conditions and Stability:	

AB-NK276-2 MARK3-BCT Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications

Product Use:	Western blotting Antibody microarrays
Antibody Dilution Recommended:	2 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals

This product is for in vitro research use only and is not intended for use in humans or animals.