

AB-NN060-7

HSP70 Antibody

Pan-specific monoclonal antibody (1.86) for monitoring the expression of bovine heat shock protein HSP70



KINEXUS

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Target Protein

Name Long:	Heat shock 70 kDa protein 1A
Alias:	HSP70 1, HSP70 2, HSP70.1, HSP72, HSPA1, HSPA1A, HSPA1B
UniProt ID:	Q27975 - Bovine
Human Predicted Mass (kDa):	70.109 (642 AA; A0A0G2JIW1-1)
Observed SDS-PAGE Mass (kDa):	70

Immunogen

Antibody Immunogen Source:	Bovine HSP70
-----------------------------------	--------------

Production

Antibody Host Species:	Mouse
Antibody Type:	Monoclonal
Antibody Ig Isotype Clone Lot:	128 IgG1
Antibody Purification:	Protein G purified
Amount:	50 µg
Antibody Concentration:	1 mg/ml
Lot Number:	15DE1
Storage Buffer:	Phosphate buffered saline, pH 7.4, 50% glycerol, 0.09% sodium azide
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles. For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

Applications

Product Use:	WB ELISA
Antibody Dilution Recommended:	WB (1:1000); optimal dilutions for assays should be determined by the user.
Antibody Species Reactivity:	Human Mouse Rat Bovine Pig
Antibody Positive Control:	1 µg/ml of SMC-113 was sufficient for detection of HSP70 in 20 µg of Hela lysate by colorimetric immunoblot analysis using goat anti-mouse IgG:HRP as the secondary antibody.
Target Detection Immunoblotting:	Detects a ~70 kDa protein.
Antibody Cross Reactivities:	Does not cross react with HSC70.

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-888-546-3987