AB-NN060-8 HSP70 Antibody

Pan-specific monoclonal antibody (5A5) for monitoring the expression of human heat shock protein HSP70



Email: info@kinexus.ca

Phone: 604-323-2547

Target Protein	Filone. 004-323-2347
Name Long:	Heat shock 70 kDa protein 1A
Alias:	HSP70 1, HSP70 2, HSP70.1, HSP72, HSPA1, HSPA1A, HSPA1B
UniProt ID:	P0DMV8 - Human
Human Predicted Mass (KDa):	70.052 (641 AA; P0DMV8-1); 63.937 (586 AA; P0DMV8-2)
Observed SDS-PAGE Mass (KDa):	70
Immunogen	
Antibody Immunogen Source:	Human recombinant HSP70 overexpressed in E.coli
Production	
Antibody Host Species:	Mouse
Antibody Type:	Monoclonal
Antibody Ig Isotype Clone Lot:	129 lgG1
Antibody Purification:	Protein G purified
Amount:	25 µg
Antibody Concentration:	1 mg/ml
Lot Number:	15DE1
Storage Buffer:	Phosphate buffered saline pH 7.4 pH7.2, 50% glycerol, 0.09% sodium azide
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.
Applications	
Product Use:	WB IHC ICC/IF IP
Antibody Dilution Recommended:	WB (1:1000), ICC/IF (1:500), IP (1 μ g); optimal dilutions for assays should be determined by the user.
Antibody Species Reactivity:	Human Mouse Rat Amphibian Chicken Fish S. cerevisiae Drosophila
Antibody Positive Control:	1 μ g/ml of SMC-162 was sufficient for detection of HSP70 in 20 μ g of heat shocked HeLa cell lysate by colorimetric immunoblot analysis using Goat antimouse IgG:HRP as the secondary antibody.
Target Detection Immunoblotting:	Detects a ~70 kDa protein.
Antibody Specificity:	Very high
Antibody Cross Reactivities:	May detect HSP70, HSC70, Grp78 and HSP72.

Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-888-546-3987