## AB-NN127 FKBP52 Antibody

Target Dretain

Pan-specific monoclonal antibody (Hi52C) for monitoring the expression of human FKBP52



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Target Protein	
Name Long:	Peptidyl-prolyl cis-trans isomerase FKBP4
Alias:	FK506 binding protein 4, FKBP4, FKBP59, HBI, HSP56, p52, p59, PPIase, Rotamase, T cell FK506 binding protein
UniProt ID:	Q02790 - Human
Human Predicted Mass (KDa):	51.805 (459 AA; Q02790-1)
Observed SDS-PAGE Mass (KDa):	52
Immunogen	
Antibody Immunogen Source:	Synthetic peptide corresponding to the residues of human FKBP52
Production	
Antibody Host Species:	Mouse
Antibody Type:	Monoclonal

Monoclonal	
87 lgG	
Protein G purified	
25 μg	
1 mg/ml	
15DE1	
Phosphate buffered saline pH 7.4, 50% glycerol, 0.09% sodium azide	
For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.	
Applications	
WB   IHC   ICC/IF   IP	

 Antibody Species Reactivity:
 Human | Mouse | Rat | Dog | Hamster

 Antibody Positive Control:
 0.5 µg/ml was sufficient for detection of FKBP52 in 20 µg total protein using WB by colorimetric immunoblot analysis using Goat Anti-Mouse IgG:HRP as the secondary.

 Target Detection Immunoblotting:
 Detects a ~52 kDa protein. Heavy chain migrates close to FKBP52 on SDS-PAGE gels.

 Antibody Specificity:
 Very High

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-888-546-3987