AB-NN181-2 PDK3-3 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human protein-serine/threonine kinase PDK3 (PDHK3)



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Target Protein

| Name Long: | Pyruvate dehydrogenase kinase isoform 3; [Pyruvate dehydrogenase [lipoamide]] kinase isozyme 3, mitochondrial |
|--------------------------------|---|
| Alias: | Kinase Pyruvate dehydrogenase kinase 3; PDK3 |
| UniProt ID: | Q15120 |
| Sequence Predicted Mass (KDa): | 48.043 (415 AA; Q15120-2); 46.939 (406 AA; Q15120) |
| Observed SDS-PAGE Mass (KDa): | 42-53 |

Immunogen

| Antibody Immunogen Source: | Human PDK3 (PDHK3) sequence peptide Cat. No.: PE-01BAE99 |
|------------------------------|---|
| Antibody Immunogen Sequence: | LEYENKSPEDPQVLDC |
| Location in Target: | Corresponds to amino acid residues L92 to D106; |
| Peptide Type: | For pan-specific recognition of target expression levels. |
| Target Phosphosite: | Not phosphorylated |

Production

| Antibody Host Species: | Rabbit |
|-----------------------------------|--|
| Antibody Type: | Polyclonal |
| Antibody Ig Isotype Clone Lot: | Immunoglobulin G |
| | The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris. |
| Production Method: | |
| Antibody Amount: | 25 µg |
| Antibody Concentration: | 1 mg/ml |
| Lot Number: | 141204 |
| Storage Buffer: | Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol |
| Storage Conditions and Stability: | For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles. |
| | |

AB-NN181-2 PDK3-3 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications

| The second secon | |
|--|---|
| Product Use: | Western blotting Antibody microarrays |
| Antibody Dilution Recommended: | 2 μg/ml for immunoblotting |
| Antibody Species Reactivity: | Human, mouse, rat and many other mammals |
| Antibody Positive Controls: | Weak reactivity with recombinant human PDK3 on protein dot blots. |
| Overall Antibody Specificity: | Medium selectivity |
| Antibody Cross Reactivities: | No immunoreactivity on protein dot blots with recombinant human PDK1, PDK2 and PDK4 |

This product is for in vitro research use only and is not intended for use in humans or animals.