# AB-NN191-4 AHA1 Antibody

Pan-specific polyclonal antibody for monitoring the expression of mouse AHA1



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

## **Target Protein**

Activator of 90 kDa heat shock protein ATPase homologue 1 (AHA1)

Activator of 90 kDa heat shock protein ATPase homologue 1, Activator of heat shock 90kDa protein ATPase homologue 1, AHA 1, AHA1, AHSA 1, AHSA1, AHSA1\_HUMAN, C14orf3, HSPC322, p38

UniProt ID:

Q8BK64 - Mouse

Human Predicted Mass (KDa):

38.274 (338 AA; O95433-1); 32.340 (288 AA; O95433-2)

Observed SDS-PAGE Mass (KDa): 38-45

### **Immunogen**

Antibody Immunogen Source: Mouse Aha1

#### Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	48 Mix
Antibody Purification:	Protein A purified
Amount:	25 μΙ
Antibody Concentration:	1 mg/ml
Lot Number:	15DE1
Storage Buffer:	Phosphate buffered saline, pH 7.4, 50% glycerol, 0.09% sodium azide
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles. For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

#### **Applications**

Product Use:	WB   IHC   ICC/IF   IP   ELISA
Antibody Dilution Recommended:	WB (1:1000); optimal dilutions for assays should be determined by the user.
<b>Antibody Species Reactivity:</b>	Human   Mouse   Rat
Antibody Positive Control:	1 $\mu$ I/ml of SPC-183 was sufficient for detection of Aha1 in 10 $\mu$ g of mixed human cell lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.
Target Detection Immunoblotting:	Detects a ~38 kDa protein. Can migrate as a 45 kDa protein on SDS-PAGE gels based on immunoreactivity with recombinant Aha1.
Antibody Specificity:	Medium-high

This product is for in vitro research use only and is not intended for use in humans or animals.