AB-NN199-1

Aquaporin 1 Antibody

Pan-specific polyclonal antibody for monitoring the expression of rat Aquaporin 1



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Target Protein

Name Long:	Aquaporin 1
Alias:	AQP1, AQP 1, AQP CHIP, AQP-1, AQP1_HUMAN, Aquaporin CHIP, Aquaporin-1, Aquaporin-CHIP, Aquaporin1, Channel forming integral protein 28kDa, Channel like integral membrane protein 28 kDa, CHIP 28, CHIP28, CO, Colton blood group, Growth factor induced delayed early response protein, MGC26324, Urine water channel, Water channel protein CHIP 29, Water channel protein CHIP29, Water channel protein for red blood cells and kidney proximal tubule
UniProt ID:	P29975 - Rat
Human Predicted Mass (KDa):	28.526 (269 AA; A0A024RA31-1)
Observed SDS-PAGE Mass (KDa):	29

Immunogen

Antibody Immunogen Source:	Produced against the N-terminal peptide (sequence N-MASEFKKKLF) of rate Aquaporin 1	t
----------------------------	---	---

Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	56 Mix
Antibody Purification:	Affinity purified
Amount:	25 μg
Antibody Concentration:	1 mg/ml
Lot Number:	15DE1
Storage Buffer:	Phosphate buffered saline pH 7.4, 50% glycerol, 0.09% sodium azide
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles. For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

Applications

Product Use:	WB IHC ICC/IF
Antibody Dilution Recommended:	WB (1:2000), ICC/IF (1:400); optimal dilutions for assays should be determined by the user.
Antibody Species Reactivity:	Human Mouse Rat
Antibody Positive Control:	0.5 μg/ml of SPC-502 was sufficient for detection of aquaporin 1 in 10 μg of rat kidney tissue lysate by colorimetric immunoblot analysis using Goat anti-rabbit lgG:HRP as the secondary antibody.
Target Detection Immunoblotting:	Detects a ~28.5 kDa protein. May detect larger glycosylated ~35-50 kDa

This product is for in vitro research use only and is not intended for use in humans or animals.