AB-NN232-1

Chapsyn-110 Antibody

Pan-specific monoclonal antibody (S18-30) for monitoring the expression of rat Chapsyn-110



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Target Protein

Name Long:

Disks large homologue 2 (Chapsyn-110)

PSD93, Chapsyn-110, MGC131811, PSD-93, Discs Large homologue 2, DLG2,
Postsynaptic density protein PSD-93, DLGH2, Channel-associated protein of synapse-110

UniProt ID: Q63622 - Rat

Human Predicted Mass (KDa): 95.000 (852 AA; B7Z2T4-1)

Observed SDS-PAGE Mass (KDa): 110

Immunogen

Antibody Immunogen Source: Fusion protein amino acids 1-852 of rat Chapsyn-110

Production

Antibody Host Species:	Mouse
Antibody Type:	Monoclonal
Antibody Ig Isotype Clone Lot:	9 lgG1
Antibody Purification:	Protein G purified
Amount:	25 μg
Antibody Concentration:	1 mg/ml
Lot Number:	15DE1
Storage Buffer:	Phosphate buffered saline, pH 7.4, 50% glycerol, 0.09% sodium azide
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles. For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

Applications

Product Use:	WB IHC IP
Antibody Dilution Recommended:	WB (1:1000), IHC (1:1000), ICC/IF (1:100); optimal dilutions for assays should be determined by the user.
Antibody Species Reactivity:	Human Mouse Rat
Antibody Positive Control:	1 μ g/ml of SMC-325 was sufficient for detection of Chapsyn -110 in 10 μ g of rat brain lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.
Target Detection Immunoblotting:	Detects a ~110 kDa protein. Weak human detection.
Antibody Specificity:	Very High
Antibody Cross Reactivities:	No cross-reactivity against other MAGUK family members, PSD95, SAP97, SAP102, expressed in transfected cells.

This product is for in vitro research use only and is not intended for use in humans or animals.