AB-NN235-1 Cpn10 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human Cpn10



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Target Protein	Canada vor 013 Phone. 604-525-2547
Name Long:	10 kDa heat shock protein, mitochondrial (cpn10)
Alias:	10kDa Chaperonin, Chaperonin 10, Cpn 10, EPF, GROES, Heat Shock 10kD protein1, HSP10, HSPE1
UniProt ID:	P61604 - Human
Human Predicted Mass (KDa):	10.932 (102 AA; P61604-1)
Observed SDS-PAGE Mass (KDa):	10
Immunogen	
Antibody Immunogen Source:	Human Cpn10 peptide AA 91-101
Production	
Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	74 Mix
Antibody Purification:	Protein A purified
Amount:	25 µg
Antibody Concentration:	1 mg/ml
Lot Number:	15DE1
Storage Buffer:	Phosphate buffered saline pH 7.4 pH 7.3, 0.02% sodium azide
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.
Applications	
Product Use:	WB IHC ICC/IF IP ELISA
Antibody Dilution Recommended:	WB (1:1000); optimal dilutions for assays should be determined by the user.
Antibody Species Reactivity:	Human Mouse Rat Bovine Dog Guinea Pig Pig Rabbit Sheep Xenopus
Antibody Positive Control:	1 μ g/ml of SPC-193 was sufficient for detection of Cpn10 in 10 μ g of Hela lysates by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.
Target Detection Immunoblotting:	Detects a ~ 10 kDa protein.
Antibody Specificity:	Very high

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-888-546-3987