AB-NN295-1

Neuroligin 3 Antibody Pan-specific monoclonal antibody (S110-29) for monitoring the

expression of rat Neuroligin 3

Torgot Drotoin



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

| Target Protein | |
|-----------------------------------|--|
| Name Long: | Neuroligin-3 |
| Alias: | Gliotactin homologue, Neuroligin-3, Nlgn3, NLGN3_HUMAN |
| UniProt ID: | Q62889 - Rat |
| Human Predicted Mass (KDa): | 93.895 (848 AA; X5DNV3-1) |
| Observed SDS-PAGE Mass (KDa): | 110 |
| | |
| Immunogen | |
| Antibody Immunogen Source: | Fusion protein amino acids 730-848 (intracellular C-terminus) of rat Neuroligin-3. Mouse: 99% identity (118/119 amino acids identical). Human: 98% identity (116/119 amino acids identical) ~60% identity with Neuroligin- 1. ~40% identity with Neuroligin-2 |
| Production | |
| Antibody Host Species: | Mouse |
| Antibody Type: | Monoclonal |
| Antibody Ig Isotype Clone Lot: | lgG1 |
| Antibody Purification: | Protein G purified |
| Amount: | 25 µg |
| Antibody Concentration: | 1 mg/ml |
| Lot Number: | 15DE1 |
| Storage Buffer: | Phosphate buffered saline pH 7.4, 50% glycerol, 0.1% sodium azide |
| Storage Conditions and Stability: | For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles. |
| Applications | |
| Product Use: | WB IHC ICC/IF |
| Antibody Dilution Recommended: | WB (1:1000); optimal dilutions for assays should be determined by the user. |
| Antibody Species Reactivity: | Human Mouse Rat |
| Antibody Positive Control: | 1 μ g/ml of SMC-471 was sufficient for detection of Neuroligin 3 in 20 μ g of rat brain lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody. |
| Target Detection Immunoblotting: | In mouse brain lysates, this antibody detects a ~110 kDa protein. |
| Antibody Specificity: | High |
| Antibody Cross Reactivities: | Does not cross-react with Neuroligin-1, -2, -4 or -4. Weak~120 KDa target protein and one very strong ~75 KDa cross-reactive protein in mouse brain lysates. |

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-888-546-3987