# AB-NN345-1

## Synaptotagmin-12 Antibody

Pan-specific monoclonal antibody (S277-7) for monitoring the expression of mouse Synaptotagmin-12



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia,

Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

#### **Target Protein**

Name Long: Synaptotagmin-12 SYT-12, SYT12, Synaptotagmin 12, Synaptotagmin XII, SytXII, Syt XII, Alias: Synaptotagmin12, SynaptotagminXII **UniProt ID:** Q920N7 - Mouse **Human Predicted Mass (KDa):** 46.537 (421 AA; Q8IV01-1)

Observed SDS-PAGE Mass (KDa):

### Immunogen

Fusion protein amino acids 168-255 (Cytoplasmic C2A domain) of mouse **Antibody Immunogen Source:** 

Synaptotagmin-12

#### **Production**

Antibody Host Species:	Mouse
Antibody Type:	Monoclonal
Antibody Ig Isotype Clone Lot:	IgG1
Antibody Purification:	Protein G purified
Amount:	25 μg
Antibody Concentration:	1 mg/ml
Lot Number:	15DE1
Storage Buffer:	Phosphate buffered saline pH 7.4, 50% glycerol, 0.09% sodium azide
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles. For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

#### **Applications**

Product Use:	WB   IHC   ICC/IF
Antibody Dilution Recommended:	WB (1:1000); optimal dilutions for assays should be determined by the user.
Antibody Species Reactivity:	Mouse   Rat
Antibody Positive Control:	1 $\mu$ g/ml of SMC-437 was sufficient for detection of Synaptotagmin-12 in 20 $\mu$ g of transiently overexpressing synaptotagmin-12 COS cell lysate by colorimetric immunoblot analysis using Goat anti-mouse IgG:HRP as the secondary antibody.
Target Detection Immunoblotting:	In mouse brain lysates, this antibody detects a ~45 kDa protein.
Antibody Specificity:	Very High
Antibody Cross Reactivities:	Does not cross-react with Synaptotagmin-6 or other Synaptotagmins.

This product is for in vitro research use only and is not intended for use in humans or animals.