

AB-NN349-1

TRAP1 Antibody

Pan-specific monoclonal antibody (3H4-2H6) for monitoring the expression of human Trap1



KINEXUS

Address: 8755 Ash Street, Suite 1
Vancouver, British Columbia,
Canada V6P 6T3

Email: info@kinexus.ca
Phone: 604-323-2547

Target Protein

Name Long:	Heat shock protein 75 kDa, mitochondrial (Trap1)
Alias:	Heat shock protein 75, Heat shock protein 75 kDa, Heat shock protein 75 kDa mitochondrial, HSP 75, HSP 90L, HSP75, HSP90L, mitochondrial, TNF receptor associated protein 1, TNFR associated protein 1, TNFR-associated protein 1, TRAP 1, TRAP-1, Trap1, TRAP1_HUMAN, Tumor necrosis factor receptor associated protein, Tumor necrosis factor type 1 receptor associated protein, Tumor necrosis factor type 1 receptor-associated protein
UniProt ID:	Q12931 - Human
Human Predicted Mass (KDa):	80.110 (704 AA; Q12931-1); 74.267 (651 AA; Q12931-2)
Observed SDS-PAGE Mass (KDa):	75

Immunogen

Antibody Immunogen Source:	Purified recombinant TRAP1
-----------------------------------	----------------------------

Production

Antibody Host Species:	Mouse
Antibody Type:	Monoclonal
Antibody Ig Isotype Clone Lot:	195 IgG1Kappa
Antibody Purification:	Protein G purified
Amount:	25 µg
Antibody Concentration:	1 mg/ml
Lot Number:	15DE1
Storage Buffer:	Phosphate buffered saline, pH 7.4, 50% glycerol, 0.09% sodium azide
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles. For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

Applications

Product Use:	WB
Antibody Dilution Recommended:	WB (1:1000); optimal dilutions for assays should be determined by the user.
Antibody Species Reactivity:	Human Rat
Antibody Positive Control:	1 µg/ml of SMC-207 was sufficient for detection of Trap-1/HSP75 in 20 µg of Human A431 lysate by ECL immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.
Target Detection Immunoblotting:	Detects a ~75 kDa protein.

This product is for in vitro research use only and is not intended for use in humans or animals.

For more information on our products please visit www.kinexusproducts.ca or contact us at 1-888-546-3987