AB-NN370-1

Synaptophysin Antibody

Pan-specific polyclonal antibody for monitoring the expression of Synaptophysin.



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Target Protein

Name Long:	Synaptophysin
Alias:	Major synaptic vesicle protein p38, SYP
UniProt ID:	P08247 - Human
Human Predicted Mass (KDa):	33.845 (313 AA; P08247-1); 20.757 (195 AA; P08247-2)

Immunogen

Antibody Immunogen Source:

Production

Antibody Host Species:	Mouse
Antibody Type:	Monoclonal
Antibody Ig Isotype Clone Lot:	Mix
Amount:	25 μg
Antibody Concentration:	1 mg/ml
Lot Number:	15DE1
Storage Buffer:	Phosphate buffered saline pH7.4, lyophilized
	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept
Storage Conditions and Stability:	at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

Applications

Product Use:	WB
Antibody Dilution Recommended:	WB (1:1000)
Antibody Species Reactivity:	Human
Target Detection Immunoblotting:	In mouse brain lysates, this antibody detects a ~38 kDa protein.
Antibody Cross Reactivities:	No other prominent cross-reactive proteins in mouse brain lysates.

This product is for in vitro research use only and is not intended for use in humans or animals.