AB-NP002-2 CDC25B-3 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human CDC25B



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Target Protein

Name Long: Cell division cycle 25B phosphatase

CDC25HU2; CDC25M2; Cell division cycle 25 B; M-phase inducer phosphatase Alias:

2; MPIP2

UniProt ID: P30305

64.987 (580 AA; P30305-1); 63.443 (566 AA; P30305-2); 60.756 (539 AA; Sequence Predicted Mass (KDa):

P30305-3); 67.241 (601 AA; P30305-4)

Observed SDS-PAGE Mass (KDa): 62-68

Immunogen

Antibody Immunogen Source: Human CDC25B sequence peptide Cat. No.: PE-01ASY75

Antibody Immunogen Sequence: CAFKDELKTFRLKTR

Location in Target: Corresponds to amino acid residues A549 to R562; Peptide Type: For pan-specific recognition of target expression levels.

Target Phosphosite: Not phosphorylated

Production

Production Method:

Antibody Host Species: Rabbit **Antibody Type:** Polyclonal

Antibody Ig Isotype Clone Lot: Immunoglobulin G

> The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by

mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The

sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH

7.0 with saturated Tris.

Antibody Amount: 25 µg

Antibody Concentration: 1 mg/ml 140724 Lot Number:

Storage Buffer: Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol

For long term storage, keep frozen at -40°C or lower. Stock solution can be kept Storage Conditions and Stability:

at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

AB-NP002-2 CDC25B-3 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

Applications

P. P. Sandara	
Product Use:	Western blotting Antibody microarrays
Antibody Dilution Recommended:	1 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Antibody Positive Controls:	Strong immunoreactivity with recombinant human CDC25B on protein dot blots.
Overall Antibody Specificity:	Medium selectivity
Antibody Cross Reactivities:	Strong immunoreactivity on protein dot blots with recombinant human CDC25C, but no immunoreactivity with recombinant human CDC25A.

This product is for in vitro research use only and is not intended for use in humans or animals.