# AB-NP039-3 DUSP5-3 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human DUSP5 (VH3)



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

## **Target Protein**

Name Long:	Dual specificity protein phosphatase 5
Alias:	Dual specificity protein phosphatase hVH3; DUS5; DUS5_HUMAN; HVH3; VH3
UniProt ID:	Q16690
Sequence Predicted Mass (KDa):	42.047 (384 AA; Q16690)
Observed SDS-PAGE Mass (KDa):	40-42

### Immunogen

Antibody Immunogen Source:	Human DUSP5 (VH3) sequence peptide Cat. No.: PE-01BAZ90
Antibody Immunogen Sequence:	CLKEAFDYIKQRRSMV
Location in Target:	Corresponds to amino acid residues L286 to V300
Peptide Type:	For pan-specific recognition of target expression levels.
Target Phosphosite:	Not phosphorylated

#### Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by
Production Method:	mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
Antibody Amount:	25 μg
Antibody Concentration:	0.29 mg/ml
Lot Number:	140724
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

## AB-NP039-3 DUSP5-3 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

## Applications

P. P	
Product Use:	Western blotting   Antibody microarrays
<b>Antibody Dilution Recommended:</b>	2 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Overall Antibody Specificity:	Medium selectivity

This product is for in vitro research use only and is not intended for use in humans or animals.