# AB-NP041-2 DUSP7-2 Antibody

Pan-specific polyclonal antibody for monitoring the expression of human DUSP7



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

## **Target Protein**

Name Long:	Dual specificity protein phosphatase 7
Alias:	Dual specificity protein phosphatase PYST2; DUS7; DUS7_HUMAN; PYST2
UniProt ID:	Q16829
Sequence Predicted Mass (KDa):	44.957 (419 AA; Q16828-1); 40.551 (368 AA; Q16828-2)
Observed SDS-PAGE Mass (KDa):	37-43

### Immunogen

Antibody Immunogen Source:	Human DUSP7 sequence peptide Cat. No.: PE-01BBE95
Antibody Immunogen Sequence:	CDNHASSEQLYFSTP
Location in Target:	Corresponds to amino acid residues C391 to P405; Pre-C-terminus
Peptide Type:	For pan-specific recognition of target expression levels.
Target Phosphosite:	Not phosphorylated

#### Production

Antibody Host Species:	Rabbit
Antibody Type:	Polyclonal
Antibody Ig Isotype Clone Lot:	Immunoglobulin G
Production Method:	The immunizing peptide was produced by solid phase synthesis on a multipep peptide synthesizer and purified by reverse-phase hplc chromatography. Purity was assessed by analytical hplc and the amino acid sequence confirmed by mass spectrometry analysis. This peptide was coupled to KLH prior to immunization into rabbits. New Zealand White rabbits were subcutaneously injected with KLH-coupled immunizing peptide every 4 weeks for 4 months. The sera from each animal was applied onto an agarose column to which the immunogen peptide was thio-linked. Antibody was eluted from the column with 0.1 M glycine, pH 2.5. Subsequently, the antibody solution was neutralized to pH 7.0 with saturated Tris.
Antibody Amount:	25 µg
Antibody Concentration:	0.5 mg/ml
Lot Number:	140724
Storage Buffer:	Phosphate buffered saline (PBS) pH7.4, 0.05% Thimerasol
Storage Conditions and Stability:	For long term storage, keep frozen at -40°C or lower. Stock solution can be kept at +4°C for more than 3 months. Avoid repeated freeze-thaw cycles.

## AB-NP041-2 DUSP7-2 Antibody



Address: 8755 Ash Street, Suite 1 Vancouver, British Columbia, Canada V6P 6T3

Email: info@kinexus.ca Phone: 604-323-2547

## Applications

Product Use:	Western blotting   Antibody microarrays
<b>Antibody Dilution Recommended:</b>	2 μg/ml for immunoblotting
Antibody Species Reactivity:	Human, mouse, rat and many other mammals
Overall Antibody Specificity:	Very high selectivity

This product is for in vitro research use only and is not intended for use in humans or animals.